APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ <b>HITACHI 911/912</b>	
TEST:	<b>PH</b>
APP. CODE:	<b>344</b>
WAVELENGTH (Sec/Pri):	<b>700 - 340</b>
ASSAY:	<b>1-POINT</b> <i>TIME: 10</i> <i>POINT: 31</i>
SAMPLE VOL:	NORMAL: <b>3</b> DECREASE: <b>2</b> INCREASE: <b>5</b>
	R1 VOLUME: <b>300</b> R2 VOLUME: <b>0</b> R3 VOLUME: <b>0</b> R4 VOLUME: <b>0</b>
ABS LIMIT:	<b>32000 - INC</b>
PROZONE LIMIT:	<b>0 - UPPER</b>
CALIB METHOD:	<b>LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)</b>
SD LIMIT:	<b>0.250</b>
DUPLICATE LIMIT:	<b>3%</b>
ST. 1 CONC:	<b>0.00</b>
EXPECTED VALUE:	<b>2.0 - 7.0</b>
UNIT:	<b>mg/dl</b>
INSTR. FACTOR (y=ax+b):	a=1    b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ <b>OLYMPUS AU 400/480/600/640/680/2700 (Test code 874)</b>	
TEST NAME:	<b>PH</b>
SAMPLE:	Volume <b>3 µl</b> Dilution <b>0 µl</b>
REAGENTS:	R1 Volume <b>300 µl</b> Dilution <b>0 µl</b> R2 Volume <b>0 µl</b> Dilution <b>0 µl</b>
WAVELENGTH:	Pri. <b>340</b> Sec. <b>700</b>
METHOD:	<b>END</b>
REACTION SLOPE:	<b>+</b>
MEASURING POINT 1:	First <b>0</b> Last <b>17</b>
MEASURING POINT 2:	First      Last
REAGENT OD LIMIT:	First L <b>-0.1</b> First H <b>0.5</b> Last L <b>-0.1</b> Last H <b>0.5</b>
DYNAMIC RANGE:	L <b>0.4</b> H <b>20</b>
CORRELATION FACTOR:	A <b>1</b> B <b>0</b>
UNIT:	<b>mg/dl</b>
CALIBRATION TYPE:	<b>AB</b>
FORMULA:	<b>Y = AX + B</b>

 Chema Diagnostica  
Via Campania 2/4  
60030 Monsano (AN) - ITALY - EU  
phone +39 0731 605064  
fax +39 0731 605672  
e-mail: mail@chema.com  
website: http://www.chema.com

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**INTENDED USE**  
Reagent for quantitative in vitro determination of phosphorus in biological fluids.

**PRINCIPLE OF THE METHOD**  
The phosphate ions react with ammonium molybdate to form a phosphomolybdate complex. The colourless phosphomolybdate complex can be measured directly by ultraviolet (UV) absorption at 340 nm. An acid pH is necessary for the formation of complexes.

**KIT COMPONENTS**

**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources.

**PHOS R1 2H500: 10 x 50 ml (liquid) white cap**  
**6U448: 8 x 56 ml (liquid) white cap**

Composizione: ammonio molibdato 0.4 mmol/l, sulphuric acid 0.21 mol/l, surfactant.

Store all components at 2-8°C.

**REAGENT PREPARATION**

Use separate reagent ready to use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.  
Caution: keep well refrigerated.

**PRECAUTIONS**

**PHOS R1:** It is not classified as hazardous.

**Standard:** It is not classified as hazardous.

**Standard:** It is not classified as hazardous.

**Standard:** It is not classified as hazardous.

**Standard:** It is not classified as hazardous.

used.  
Inorganic phosphate in whole blood specimens may either decrease or increase with time, depending on the type of specimen, temperature, and duration of storage. Levels in plasma or serum are increased by prolonged storage with cells at room temperature or 37 °C; it is important to promptly separate serum or plasma from erythrocytes. Hemolyzed specimens are unacceptable because erythrocytes contain high concentrations of organic phosphate esters, which can be hydrolyzed to inorganic phosphate during storage. Inorganic phosphate increases by 4 to 5 mg/dl per day in hemolyzed specimens stored at 4°C. Glucose phosphate, creatine phosphate, and other organic phosphates may also be hydrolyzed by assay conditions, resulting in overestimation of inorganic phosphate levels. Phosphate is considered to be stable in serum that has been separated from the clot for days at 4°C and months when frozen.

Urine samples should be collected in 6 mol/l HCl, 20-30 ml for a 24 hours specimen, to avoid precipitation of phosphate complexes.  
Dilute urine samples 1:20 with purified water before assay.

**EXPECTED VALUES**  
serum/plasma (adults): 2.5 - 4.5 mg/dl (0.81 - 1.45 mmol/l)  
serum/plasma (children): 4.0 - 7.0 mg/dl (1.29 - 2.26 mmol/l)

urine (nonrestricted diet): 0.4 - 1.3 g/24h (12.9 - 42.2 mmol/24h)

Each laboratory should establish appropriate reference intervals related to its population.

**QUALITY CONTROL AND CALIBRATION**

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

**QUANTINORM CHEMA - MULTINORM CHEMA** with normal or close to normal control values  
**QUANTIPATH CHEMA - MULTIPATH CHEMA** with pathological control values.

If required, a multiparametric, human based calibrator is available:

**AUTOCAL H**

Please contact Customer Care for further information.

**TEST PERFORMANCE**

**Linearity**  
the method is linear up to 20 mg/dl.  
If the limit value is exceeded, it is suggested to dilute sample 1+9 with distilled water and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**  
the limit of detection is 0.4 mg/dl.

**Interferences**  
no interference was observed by the presence of:  
bilirubin ≤ 25 mg/dl  
hemoglobin ≤ 100 mg/dl  
Hemolysis interferes.  
Both positive and negative interferences with lipemic samples has been observed.

**Precision**  
intra-assay (n=10)    mean (mg/dl)    SD (mg/dl)    CV%  
sample 1      3.50      0.05      1.50  
sample 2      5.87      0.11      1.90

inter-assay (n=20)    mean (mg/dl)    SD (mg/dl)    CV%  
sample 1      3.41      0.08      2.40  
sample 2      5.84      0.11      1.90

**Methods comparison**  
a comparison between Chema and a commercially available product gave the following results:

phosphorus Chema = x  
phosphorus competitor = y  
n = 102  
  
y = 1.005x - 0.109 mg/dl      r<sup>2</sup> = 0.975

**WASTE DISPOSAL**

This product is made to be used in professional laboratories.

P501: Dispose of contents according to national/international regulations.