APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912		
TEST:	LDL	
APP. CODE:	059	
WAVELENGTH (Sec/Pri):	700 - 600	
ASSAY:	2-POINT END	TIME: 10 POINT: 16 - 31
SAMPLE VOL:	NORMAL: 3 DECREASE: 2 INCREASE: 5	
	R1 VOLUME: 270 R2 VOLUME: 0 R3 VOLUME: 90 R4 VOLUME: 0	
ABS LIMIT:	32000 - INC	
PROZONE LIMIT:	0 - UPPER	
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)	
SD LIMIT:	0.250	
DUPLICATE LIMIT:	3%	
ST. 1 CONC:	0.0	
EXPECTED VALUE:	70.0 - 220.0	
UNIT:	mg/dl	
INSTR. FACTOR (y=ax+b):	a=1 b= 0	

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 867)		
TEST NAME:	LDL	
SAMPLE:	Volume 3 µl	Dilution 0 µl
REAGENTS:	R1 Volume 270 µl R2 Volume 90 µl	Dilution 0 µl Dilution 0 µl
WAVELENGHT:	Pri. 600 Sec. 700	
METHOD:	END	
REACTION SLOPE:	+	
MEASURING POINT 1:	First 0 Last 27	
MEASURING POINT 2:	First 0 Last 10	
REAGENT OD LIMIT:	First L -0.1 First H 0.5 Last L -0.1 Last H 0.5	
DYNAMIC RANGE:	L 1.0 H 500	
CORRELATION FACTOR:	A 1 B 0	
UNIT:	mg/dl	
CALIBRATION TYPE:	AB	
FORMULA:	Y = AX + B	

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

ITALIANO rev. 09/05/2019

LDL-direct FL		
DL 2H080	3 x 20 + 1 x 20 ml	
DL 6U120	2 x 45 + 2 x 15 ml	

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INTENDED USE
Reagent for quantitative in vitro determination of LDL-cholesterol in biological fluids.

PRINCIPLE OF THE METHOD
When a sample is mixed with Reagent R1, the protecting reagent binds to LDL and protects LDL from enzyme reactions. Cholesterol esterase (CHE) and cholesterol oxidase (CO) react with non-LDL lipoproteins [chylomicrons (CM), very low density lipoproteins (VLDL) and HDL]. Hydrogen peroxide produced by the enzyme reactions with non-LDL cholesterol is decomposed by catalase in Reagent R1. When Reagent R2 is added, the protecting reagent is removed from LDL and catalase is inactivated. In this second process, CHE and CO react only with LDL-C. Hydrogen peroxide produced by the enzyme reactions with LDL-C yields a color complex upon oxidase condensation with N-(3-sulfopropyl)-3-methoxy-5-methylaniline (HMMPs) and 4-aminoantipyrine (4AA) in the presence of peroxidase (POD). By measuring the absorbance of the blue color complex produced at approximately 600 nm, the LDL cholesterol concentration in the sample can be calculated.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label. Keep away from direct light sources.

LDL-C R1 2H080: 3 x 20 ml (liquid) white cap
6U120: 2 x 45 ml (liquid) white cap

Composition: Good's buffer pH 7.0, cholesterol esterase, cholesterol oxidase, HMMPs and catalase.

LDL-C R2 2H080: 1 x 20 ml (liquid) red cap
6U120: 2 x 15 ml (liquid) red cap

Composition: Good's buffer pH 7.0, 4-aminoantipyrine, POD.

Store at 2-8°C. Do not freeze.

REAGENT PREPARATION

Use separate reagents ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: 30 days at 2-8°C.

PRECAUTIONS
LDL-C R1: **Warning.** May cause an allergic skin reaction (H317). Wear protective gloves/protective clothing/eye protection/face protection (P280). IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower] (P303+P361+P353). If skin irritation or rash occurs: Get medical advice (P333+P313).

LDL-C R2: It is not classified as hazardous.

N-acetylcysteine (NAC), metamazole and acetaminophen may cause interference in the Trinder reaction.^(4,5)
To avoid interference, the blood withdrawal should be performed before drug administration.

SPECIMEN
Serum, plasma heparinate.
Anticoagulants such as heparin, citrate, oxalate and EDTA do not have significant influences on the assay when they are used in their usual amounts.
Samples with triglyceride concentrations exceeding 1000 mg/dl should be diluted an reanalyzed

Use fresh specimens. Do not use specimens that repeated freeze-thaw, because lipoproteins may be denatured.

EXPECTED VALUES
normal values: 76 - 218 mg/dl
NCEP ATP's Decision cut-off points for LDL-C:
desirable: < 130 mg/dl
borderline high risk for CHD 130 - 159 mg/dl
high risk for CHD ≥ 160 mg/dl

QUALITY CONTROL AND CALIBRATION
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

QUANTINORM CHEMA
with normal or close to normal control values
QUANTIPATH CHEMA
with pathological control values.
If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE
Linearity
the method is linear up to 500 mg/dl.
If the limit is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
the limit of detection is 1 mg/dl.

Interferences
no interference was observed by the presence of:
hemoglobin ≤ 500 mg/dl
bilirubin (free) ≤ 50 mg/dl
bilirubin (conjugated) ≤ 40 mg/dl
ascorbic acid ≤ 50 mg/dl
lipid ≤ 2000 mg/dl

Precision
intra assay (n=21) mean (mg/dl) SD (mg/dl) CV%
sample 1 85.6 1.46 1.71
sample 2 129.6 2.28 1.76
inter-assay (n=9) mean (mg/dl) SD (mg/dl) CV%
sample 3 87.8 1.69 1.92
sample 4 129.6 1.99 1.53

Methods comparison
a comparison between LDL-direct FL and a CDC Reference method (beta-quantification) gave the following results:

LDL-direct FL Chema = x
Reference method = y
n = 25
y = 1.0015x - 0.715 mg/dl r² = 0.986

WASTE DISPOSAL
This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.

LDL-direct FL		
DL 2H080	3 x 20	+ 1 x 20 ml
DL 6U120	2 x 45	+ 2 x 15 ml

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**BIBLIOGRAFIA / BIBLIOGRAPHY /
 BIBLIOGRAPHIE / BIBLIOGRAFÍA /
 БИБЛИОГРАФИЯ**

1) Burtis, C. A and Ashwood, E. R., Ed. Tietz Textbook of
 Clinical Chemistry, 2nd Ed., Saunders, Philadelphia, 1994.
 2) NIH Publication No 95-3044, Recommendations on
 Lipoprotein Measurement (1995).
 3) Japan Atherosclerosis Society: Guidelines for Pre-
 vention of Atherosclerotic Cardiovascular Diseases, 5-7
 (2002).
 4) N-acetylcysteine interference of Trinder-based assays.
 Genzen JR, Hunsaker JJ, Nelson LS, Faine BA, Krasowski
 MD. Clin Biochem. 2016 Jan;49(1-2):100-4
 5) Drug interference in Trinder reaction.
 Wiewiorka O, Čermáková Z, Dastych M. Euromedlab
 2017. ISSN 1437-4431

IVD	dispositivo medico-diagnostico <i>in vitro</i> <i>in vitro</i> diagnostic medical device dispositif médical de diagnostic <i>in vitro</i> producto sanitario para diagnóstico <i>in vitro</i> <i>in vitro</i> diagnostischeeske medischeeske toestel
LOT	numero di lotto batch code numéro de lot número de lote лот выпуска
REF	numero di catalogo catalogue number référence catalogue número de catálogo номер по каталогу
	limite di temperatura temperature limit limite de température limite de temperatura диапазон температуры при хранении
	usare entro la data use-by date utiliser avant la date utilizar por fecha срок годности
	attenzione caution attention atención внимание
	consultare le istruzioni d'uso consult instructions for use consulter les instructions d'utilisation consultar las instrucciones de uso смотреть рабочие инструкции