APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	<b>RF</b>
APP. CODE:	<b>358</b>
WAVELENGTH (Sec/Pri):	<b>/ - 340</b>
ASSAY:	<b>2 POINT END</b> <i>TIME: 10</i> <i>POINT: 16-31</i> <i>DILUENT: Water</i>
SAMPLE VOL:	NORMAL: <b>16</b> DECREASE: <b>8</b> INCREASE: <b>25</b>
	R1 VOLUME: <b>250</b> R2 VOLUME: <b>0</b> R3 VOLUME: <b>50</b> R4 VOLUME: <b>0</b>
ABS LIMIT:	<b>32000 - INC</b>
PROZONE LIMIT:	<b>0 - UPPER</b>
CALIB METHOD:	<b>LOGIT-LOG (5P)</b> <i>(POINT: 5 - SPAN: 5 - WEIGHT: 0)</i>
SD LIMIT:	<b>0.250</b>
DUPLICATE LIMIT:	<b>3%</b>
ST. 1 CONC:	<b>0.0</b>
EXPECTED VALUE:	<b>3 - 20</b>
UNIT:	<b>IU/ml</b>
INSTR. FACTOR (y=ax+b):	a=1    b=0

**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ  
OLYMPUS AU 400/480/600/640/680/2700 (Test code 891)**

TEST NAME:	<b>RF</b>	
SAMPLE:	Volume <b>16 µl</b>	Dilution <b>0 µl</b>
REAGENTS:	R1 Volume <b>250 µl</b> R2 Volume <b>50 µl</b>	Dilution <b>0 µl</b> Dilution <b>0 µl</b>
WAVELENGTH:	Pri. <b>340</b> Sec. -	
METHOD:	<b>END</b>	
REACTION SLOPE:	<b>+</b>	
MEASURING POINT 1:	First <b>0</b>	Last <b>27</b>
MEASURING POINT 2:	First <b>0</b>	Last <b>10</b>
REAGENT OD LIMIT:	First L <b>-0.1</b> Last L <b>-0.1</b>	First H <b>0.5</b> Last H <b>0.5</b>
DYNAMIC RANGE:	L <b>3</b>	H <b>500</b>
CORRELATION FACTOR:	A <b>1</b>	B <b>0</b>
UNIT:	<b>IU/ml</b>	
CALIBRATION TYPE:	<b>5AB</b>	
FORMULA:	<b>POLYGONAL</b>	

Le modifiche sono indicate con una linea verticale sul lato del paragrafo interessato.  
Addition, deletions or changes are indicated with a vertical line on the side of the affected paragraph.  
Les ajouts, suppressions ou modifications par une ligne verticale sur le côté du paragraphe concerné.  
Las adiciones, eliminaciones o cambios se indican con una línea vertical al lado del párrafo correspondiente.  
Изменения обозначаются вертикальной линией сбоку от соответствующего абзаца.

 Chema Diagnostica Srl  
Via Campania 2/4  
60030 Monsano (AN) - ITALY - EU  
phone +39 0731 605064  
fax +39 0731 605672  
e-mail: mail@chema.com  
website: http://www.chema.com

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Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

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**INTENDED USE**

Reagent for quantitative in vitro determination of rheumatoid factor in biological fluids (serum) and intended as aid to diagnosis and prognosis of Rheumatoid Arthritis and autoimmune diseases<sup>1</sup>.  
The IVD devie can be used both in manual or on automatic analyzers. The product is intended for professional use in clinical laboratories.

**TEST PRINCIPLE**

Rheumatoid Factor (RF) reacts, via antigen-antibody reaction, with aggregated human IgG. The turbidity produced in this reaction is proportional to the concentration of RF in the sample, and it can be measured at the wavelength of 340 nm<sup>2-4</sup>.

**MATERIALS PROVIDED AND COMPOSITION**

**RF R1**      2H100: 2 x 45 ml (liquid) white cap  
              6U140: 2 x 45 ml (liquid) white cap

Composition: Good's buffer, stabilizers and preservatives.

**RF R2**      2H100: 2 x 9 ml (liquid) red cap  
              6U140: 2 x 9 ml (liquid) red cap

Composition: Heat-aggregated human IgG ≤ 0.5 mg/ml, stabilizers and preservatives.

**MATERIALS REQUIRED BUT NOT SUPPLIED**

Current laboratory instrumentation. Automatic micro-pipettes. Glass or high quality polystyrene cuvettes. Saline solution.

**REAGENT PREPARATION**

Use ready-to-use reagents.

**STABILITY AND STORAGE**

Store all kit components at 2-8°C.  
Stability of reagents: up to expiration date claimed on label storing at 2-8°C;  
Stability after first opening of reagent bottle: preferably within 60 days storing at 2-8°C;  
On-board stability on Hitachi: up to 4 weeks;  
On-board stability on Olympus: up to 4 weeks.

**PRECAUTIONS**

**RF R1:** It is not classified as hazardous.

**RF R2:** It is not classified as hazardous.

Follow required safety procedures when handling all laboratory reagents.

**SPECIMEN**

Serum.  
Samples are stable for 8 days when stored at 2-8°C and for 3 months at -20°C<sup>5</sup>.

**REFERENCE INTERVALS**

Adults<sup>1</sup> < 20 IU/ml

Each laboratory should establish its own reference ranges based on its population.

**QUALITY CONTROL AND CALIBRATION**

Calibration is required with each change in reagent lot number. It is suggested to verify calibration with at least one level of an internal quality control. If control results fall outside acceptable ranges, recalibration may be necessary. For this purpose the following human based control sera are available:

**RHEUMATOID FACTOR CONTROL SET**

For calibration, use the following product:

**RF STANDARD**

Please contact Customer Care for further information.

**TEST PERFORMANCE**

**Sensitivity / Limit of Detection (LOD)**  
The method is able to detect up to 3.0 IU/ml.

**Analytical specificity:**

**Interferences**  
interference does not occur with the following:  
Heparin: ≤ 50 mg/l  
Natrium Citraat ≤ 1000 mg/dl  
Hemoglobin ≤ 1000 mg/dl  
Bilirubin ≤ 30 mg/dl  
Intralipid ≤ 2500 mg/dl  
Ascorbic Acid ≤ 50 mg/dl  
EDTA ≤ 5 mg/d

**Veridicity<sup>6</sup>**  
BIAS% < 6.36

**Accuracy:**  
**Trueness<sup>6</sup>**  
Total observed error% < 13.50 (allowable total error)

<b>Precision</b>			
<i>Repeatability</i> (intra-assay)			
n = 20	mean (IU/ml)	SD (IU/ml)	CV%
sample 1	29.2	1.06	3.65
sample 2	105.7	2.85	2.69
sample 3	204.0	3.13	1.54

<i>Reproducibility</i> (inter-assay)			
days = 12	mean (IU/ml)	SD (IU/ml)	CV%
sample 1	24.6	0.88	3.57
sample 2	102.2	1.37	1.34
sample 3	190.2	3.63	1.91

**Measurement range**  
The lower limit is 10.0 mg/l<sup>7</sup>.