**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ  
HITACHI 911/912**

TEST:	<b>AMY</b>	
APP. CODE:	<b>334</b>	
WAVELENGTH (Sec/Pri):	<b>700 - 415</b>	
ASSAY:	<b>RATE-A</b>	<b>TIME: 10 POINT: 7 - 14</b>
SAMPLE VOL:	NORMAL: 7 DECREASE: 5 INCREASE: 10	
	R1 VOLUME: <b>280</b> R2 VOLUME: <b>0</b> R3 VOLUME: <b>0</b> R4 VOLUME: <b>0</b>	
ABS LIMIT:	<b>32000 - INC</b>	
PROZONE LIMIT:	<b>0 - UPPER</b>	
CALIB METHOD:	<b>LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)</b>	
SD LIMIT:	<b>0.250</b>	
DUPLICATE LIMIT:	<b>3%</b>	
ST. 1 CONC:	<b>0.00</b>	
EXPECTED VALUE:	<b>0 - 96</b>	
UNIT:	<b>U/l</b>	
INSTR. FACTOR (y=ax+b):	<b>a=1 b=0</b>	

**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ  
OLYMPUS AU 400/480/600/640/680/2700 (Test code 859)**

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**INTENDED USE**

Reagent for quantitative in vitro determination of amylase in biological fluids.

**PRINCIPLE OF THE METHOD**

The enzyme α-amylase hydrolyzes the 2-chloro-4-nitrophenyl-α-D-maltotrioside (CNP-G3) to release 2-chloro-4-nitrophenol and form 2-chloro-4-nitrophenyl-α-D-maltoside (CNP-G2), maltotriose (G3) and glucose (G). The rate of formation of the 2-chloro-4-nitrophenol can be detected spectrophotometrically at 405 nm to give a measurement of α-amylase activity in the sample.

**KIT COMPONENTS**

**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources.

**AMY R1 2H120: 6 x 20 ml (liquid) white cap  
6U168: 3 x 56 ml (liquid) white cap**

Composition: CNP-G3 2.3 mM, NaCl 350 mM, calcium acetate 6 mM, potassium thiocyanate 600 mM, Good's buffer pH 6.0 100 mM, stabilizers and non-reactive components.

Store all components at 2-8°C.

**REAGENT PREPARATION**

Use reagent ready to use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: preferably within 60 days at 2-8°C.

**PRECAUTIONS**

**AMY R1: Danger.** Causes serious eye damage (H318).  
 Wear protective gloves. Eye protection (P280). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). Immediately call a doctor (P310).

**SPECIMEN**

Serum, plasma (heparinate only). Urine.  
Amylase is stable in serum and plasma sample up to 2 months at 2-8°C.

**EXPECTED VALUES**

Serum - plasma: < 96 U/l (< 1.60 µkat/l)  
Urine: < 480 U/l (< 8.00 µkat/l)

Each laboratory should establish appropriate reference intervals related to its population.

**QUALITY CONTROL AND CALIBRATION**

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

**QUANTINORM CHEMA**  
with normal or close to normal control values  
**QUANTIPATH CHEMA**  
with pathological control values.  
If required, a multiparametric, human based calibrator is available:  
**AUTOCAL H**

Please contact Customer Care for further information.

**TEST PERFORMANCE**

**Linearity**  
the method is linear up to 2000 U/l.  
If a ΔA/min of 0.500 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**  
the limit of detection is 0.91 U/l.

**Interferences**  
no interference was observed by the presence of:  
hemoglobin ≤ 500 mg/dl  
bilirubin ≤ 50 mg/dl  
lipids ≤ 1200 mg/dl

**Precision**

intra-assay (n=10)	mean (U/l)	SD (U/l)	CV%
sample 1	67.89	0.97	1.42
sample 2	171.67	2.61	1.52

inter-assay (n=21)	mean (U/l)	SD (U/l)	CV%
sample 1	67.81	1.93	2.85
sample 2	175.16	4.92	2.81

**Methods comparison**  
a comparison between Chema and a commercially available product gave the following results:

$$\begin{aligned} \text{Amylase Chema} &= x \\ \text{Amylase competitor} &= y \\ n &= 155 \\ y &= 1.071x - 0.54 \text{ U/l} \quad r^2 = 0.997 \end{aligned}$$

**WASTE DISPOSAL**

This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.