APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ПРОГРАММА HITACHI 911/912	
TEST:	RF
APP. CODE:	358
WAVELENGTH (Sec/Pri):	/ - 340
ASSAY:	2 POINT END
	TIME: 10
	POINT: 16-31
	DILUENT: Water
SAMPLE VOL:	NORMAL: 16
	DECREASE: 8
	INCREASE: 25
	R1 VOLUME: 250
	R2 VOLUME: 0
	R3 VOLUME: 50
	R4 VOLUME: 0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LOGIT-LOG (5P)
	(POINT: 5 - SPAN: 5 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.0
EXPECTED VALUE:	3 - 20
UNIT:	IU/ml
INSTR. FACTOR (y=ax+b):	a=1 b= 0

**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ПРОГРАММА
OLYMPUS AU 400/480/600/640/680/2700 (Test code 891)**

TEST NAME:	RF	
SAMPLE:	Volume 16 µl	Dilution 0 µl
REAGENTS:	R1 Volume 250 µl	Dilution 0 µl
	R2 Volume 50 µl	Dilution 0 µl
WAVELENGTH:	Pri. 340 Sec. -	
METHOD:	END	
REACTION SLOPE:	+	
MEASURING POINT 1:	First 0	Last 27
MEASURING POINT 2:	First 0	Last 10
REAGENT OD LIMIT:	First L -0.1	First H 0.5
	Last L -0.1	Last H 0.5
DYNAMIC RANGE:	L 3	H 500
CORRELATION FACTOR:	A 1	B 0
UNIT:	IU/ml	
CALIBRATION TYPE:	5AB	
FORMULA:	POLYGONAL	

Le modifiche sono indicate con una linea verticale sul lato del paragrafo interessato.
Addition, deletions or changes are indicated with a vertical line on the side of the affected paragraph.
Les ajouts, suppressions ou modifications par une ligne verticale sur le côté du paragraphe concerné.
Las adiciones, eliminaciones o cambios se indican con una línea vertical al lado del párrafo correspondiente.
Изменения обозначаются вертикальной линией сбоку от соответствующего абзаца.

 Chema Diagnostica Srl
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

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Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

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INTENDED USE

Reagent for quantitative in vitro determination of rheumatoid factor in biological fluids (serum) and intended as aid to diagnosis and prognosis of Rheumatoid Arthritis and autoimmune diseases¹.
The IVD devie can be used both in manual or on automatic analyzers. The product is intended for professional use in clinical laboratories.

TEST PRINCIPLE

Rheumatoid Factor (RF) reacts, via antigen-antibody reaction, with aggregated human IgG. The turbidity produced in this reaction is proportional to the concentration of RF in the sample, and it can be measured at the wavelength of 340 nm²⁻⁴.

MATERIALS PROVIDED AND COMPOSITION

RF R1 2H100: 2 x 45 ml (liquid) white cap
6U140: 2 x 45 ml (liquid) white cap

Composition: Good's buffer, stabilizers and preservatives.

RF R2 2H100: 2 x 9 ml (liquid) red cap
6U140: 2 x 9 ml (liquid) red cap

Composition: Heat-aggregated human IgG ≤ 0.5 mg/ml, stabilizers and preservatives.

MATERIALS REQUIRED BUT NOT SUPPLIED

Current laboratory instrumentation. Automatic micro-pipettes. Glass or high quality polystyrene cuvettes. Saline solution.

REAGENT PREPARATION

Use ready-to-use reagents.

STABILITY AND STORAGE

Store all kit components at 2-8°C.
Stability of reagents: up to expiration date claimed on label storing at 2-8°C;
Stability after first opening of reagent bottle: preferably within 60 days storing at 2-8°C;
On-board stability on Hitachi: up to 4 weeks;
On-board stability on Olympus: up to 4 weeks.

PRECAUTIONS

RF R1: It is not classified as hazardous.

RF R2: It is not classified as hazardous.

Follow required safety procedures when handling all laboratory reagents.

SPECIMEN

Serum.
Samples are stable for 8 days when stored at 2-8°C and for 3 months at -20°C⁵.

REFERENCE INTERVALS

Adults¹ < 20 IU/ml

Each laboratory should establish its own reference ranges based on its population.

QUALITY CONTROL AND CALIBRATION

Calibration is required with each change in reagent lot number. It is suggested to verify calibration with at least one level of an internal quality control. If control results fall outside acceptable ranges, recalibration may be necessary. For this purpose the following human based control sera are available:

RHEUMATOID FACTOR CONTROL SET

For calibration, use the following product:

RF STANDARD

Please contact Customer Care for further information.

TEST PERFORMANCE

Sensitivity / Limit of Detection (LOD)
The method is able to detect up to 3.0 IU/ml.

Analytical specificity:

Interferences
interference does not occur with the following:
Heparin: ≤ 50 mg/l
Natrium Citraat ≤ 1000 mg/dl
Hemoglobin ≤ 1000 mg/dl
Bilirubin ≤ 30 mg/dl
Intralipid ≤ 2500 mg/dl
Ascorbic Acid ≤ 50 mg/dl
EDTA ≤ 5 mg/d

Veridicity⁶
BIAS% < 6.36

Accuracy:
Trueness⁶
Total observed error% < 13.50 (allowable total error)

Precision			
<i>Repeatability</i> (intra-assay)			
n = 20	mean (IU/ml)	SD (IU/ml)	CV%
sample 1	29.2	1.06	3.65
sample 2	105.7	2.85	2.69
sample 3	204.0	3.13	1.54

<i>Reproducibility</i> (inter-assay)			
days = 12	mean (IU/ml)	SD (IU/ml)	CV%
sample 1	24.6	0.88	3.57
sample 2	102.2	1.37	1.34
sample 3	190.2	3.63	1.91

Measurement range

The lower limit is 10.0 mg/l⁷.