

COPPER

CU 0100 CH4 x 25 ml

INTENDED USE

Reagent for quantitative in vitro determination of copper in biological fluids.

SUMMARY OF TEST

The major functions of copper metalloproteins involve oxidation-reduction reactions; most known copper-containing enzymes bind and react directly with molecular oxygen. Copper is an integral component of many metallo-enzymes, including ceruloplasmin, cytochrome c oxidase, superoxide dismutase, dopamine-β-hydroxylase, ascorbate oxidase, lysyl oxidase, and tyrosinase.

PRINCIPLE OF THE METHOD

3,5-Di-Br-PAESA combines with Cu(II) to form a blue-violet complex, the absorbance of which is measured at 580 nm. The reaction has high specificity and interference from other cations is avoided, due to specific pH and environment.

KIT COMPONENTS

For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label at 2-8°C.
Keep away from direct light sources.

CU R1 2 x 25 ml (liquid) blue cap
Composition: acetate buffer 100 mM pH 4.90, surfactants and preservatives.

CU R2 2 x 25 ml (liquid) red cap
Composition: 3,5 Di-Br-PAESA 10 mM.

Standard: copper solution 200 µg/dl - 5 ml

Store all components at 2-8°C.

MATERIALS REQUIRED BUT NOT SUPPLIED

Current laboratory instrumentation. Spectrophotometer UV/VIS with thermostatic cuvette holder. Automatic micropipettes. Glass or high quality polystyrene cuvettes. Saline solution.

REAGENT PREPARATION

Mix reagents R1 and R2 in equal parts.
Caution! Reagent R1 may precipitate during refrigerated storage. It is recommended to let it redissolve at room temperature before use. Mix thoroughly after complete dissolution.
Reagent R1 may also precipitate after mixing with R2 if the mixture is stored refrigerated. If precipitation occurs, allow the reaction solution (R1 + R2) to reach room temperature and mix thoroughly until the precipitate is completely dissolved.

Stability of working reagent: 30 days at 2-8°C and 7 days at room temperature, well closed.
Stability of unopened vials: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C.

PRECAUTIONS

Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.
Perform the test according to the general “Good Laboratory Practice” (GLP) guidelines.

SPECIMEN

Serum (preferred), plasma heparinate.
Copper is stable 7 days at 2-8°C and 1 month at -20°C.

TEST PROCEDURE			
Wavelength:	580 nm (allowed 570 ÷ 600 nm)		
Lightpath:	1 cm		
Temperature:	25, 30 or 37°C		
dispense:	blank	standard	sample
reagent	1.5 ml	1.5 ml	1.5 ml
water	100 µl	-	-
standard	-	100 µl	-
sample	-	-	100 µl
Mix, incubate at 25, 30 or 37°C for 5 minutes. Read absorbances of standard (As) and samples (Ax) against reagent blank.			

RESULTS CALCULATION

serum/plasma sample:
copper µg/dl = Ax/As x 200 (standard value)

RESULTS CALCULATION

serum/plasma sample:

copper µg/dl = Ax/As x 200 (standard value)

EXPECTED VALUES

men:	70 - 140 µg/dl	(11.0 - 22.0 µmol/l)
women:	80 - 155 µg/dl	(12.6 - 24.4 µmol/l)
pregnant women:	118 - 302 µg/dl	(18.5 - 47.4 µmol/l)
children 6-12 y:	80 - 190 µg/dl	(12.6 - 29.9 µmol/l)
infants:	20 - 70 µg/dl	(3.1 - 11.0 µmol/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

QUANTINORM CHEMA
with normal or close to normal control values

QUANTIPATH CHEMA
with pathological control values.
Please contact Customer Care for further information.

TEST PERFORMANCE

Linearity
the method is linear up to 500 µg/dl.
If the limit value is exceeded, it is suggested to dilute sample 1+9 with distilled water and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
the limit of detection is 4 µg/dl.

Interferences
no interference was observed by the presence of:
hemoglobin ≤ 120 mg/dl
bilirubin ≤ 30 mg/dl
Lipids interfere.

Precision				
intra-assay (n=10)	mean (µg/dl)	SD (µg/dl)	CV%	
sample 1	120.00	3.06	2.50	
sample 2	268.50	3.14	1.20	

	inter-assay (n=20)	mean (µg/dl)	SD (µg/dl)	CV%
sample 1		120.99	3.36	2.80
sample 2		265.19	5.73	2.22

Methods comparison
a comparison between Chema and a commercially available product gave the following results:

Copper Chema = x
Copper competitor = y
n = 82

$y = 1.046x - 6.67 \text{ µg/dl}$ $r^2 = 0.984$

WASTE DISPOSAL

This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.

REFERENCES

K.Ueno, T.Imamura, K.L.Cheng - Handbook of organic analytical reagents - CRC Press (1992).
Clin.Chem. 35/4, 552-554 (1989)
Tietz Textbook of Clinical Chemistry, Second Edition, Bur-tis-Ashwood (1994).

MANUFACTURER

Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

SYMBOLS

IVD

in vitro diagnostic medical device

LOT

batch code

REF

catalogue number

temperature limit

use-by date

caution

consult instructions for use