APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	GOT
APP. CODE:	327
WAVELENGTH (Sec/Pri):	700 - 340
ASSAY:	RATE-A <i>TIME: 10</i> <i>POINT: 20-30</i> <i>DILUENT: water</i>
SAMPLE VOL:	NORMAL: 25 DECREASE: 20 INCREASE: 30
	R1 VOLUME: 200 <i>DILUENT: 5</i> R2 VOLUME: 0 R3 VOLUME: 50 <i>DILUENT: 5</i> R4 VOLUME: 0
ABS LIMIT:	6500 - DEC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.0
EXPECTED VALUE:	0.0 - 35.0
UNIT:	U/l
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 853)	
TEST NAME:	GOT
SAMPLE:	Volume 25 µl Dilution 0 µl
REAGENTS:	R1 Volume 200 µl Dilution 0 µl R2 Volume 50 µl Dilution 0 µl
WAVELENGTH:	Pri. 340 Sec. 700
METHOD:	RATE
REACTION SLOPE:	-
MEASURING POINT 1:	First 14 Last 23
MEASURING POINT 2:	First Last
REAGENT OD LIMIT:	First L 0.8 First H 2.0 Last L 0.8 Last H 2.0
DYNAMIC RANGE:	L 0.5 H 440
CORRELATION FACTOR:	A 1 B 0
UNIT:	U/l
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

ITALIANO rev. 26/09/2016

GOT/AST FL IFCC	
GO 2H400	8 x 40 + 4 x 20 ml
GO 6U420	6 x 56 + 6 x 14 ml

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INTENDED USE
Reagent for quantitative in vitro determination of GOT in biological fluids.

PRINCIPLE OF THE METHOD
The enzyme aspartate aminotransferase (EC 2.6.1.1; L-Aspartate:2-Oxoglutarate Aminotransferase, AST or AspAT; Glutamate Oxaloacetate Transaminase, GOT) catalyzes the transaminase reaction between L-Aspartate and 2-Oxoglutarate. The 2-Oxalacetate formed, is reduced to malate in the presence of MDH. As the reactions proceed, NADH is oxidized to NAD. The disappearance of NADH per unit time is followed by measuring the decrease in absorbance at 340 nm. The present method has been made according to IFCC (2002).

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

GOT R1 **2H400: 8 x 40 ml (liquid) white cap**
6U420: 6 x 56 ml (liquid) white cap

GOT R2 **2H400: 4 x 20 ml (liquid) red cap**
6U420: 6 x 14 ml (liquid) red cap

Composition in the test: Tris buffer 80 mM pH 7.65, L-aspartate 240 mM, 2-Oxoglutarate 12 mM, NADH 0.18 mM, MDH ≥ 600 U/l, LDH ≥ 900U/l.

Store all components at 2-8°C.

REAGENT PREPARATION
Use separate reagent ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.
Caution: keep well refrigerated.

PRECAUTIONS
GOT R1: Warning. Causes serious eye irritation (H319). Causes skin irritation (H315). Wear protective gloves. Eye protection (P280). IF ON SKIN: Wash with plenty of water (P302+P352). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). If eye irritation persists: get medical advice (P337+P313).

GOT R2: It is not classified as hazardous.

SPECIMEN

Serum, plasma.
Specimens should be protected from direct exposure to light. Samples stored at 2-8°C in the dark are stable up to 3 days and 1 month at -20°C.

EXPECTED VALUES
Men: < 35 U/l (< 0.58 µkat/l)
Women: < 31 U/l (< 0.52 µkat/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

QUANTINORM CHEMA with normal or close to normal control values
QUANTIPATH CHEMA with pathological control values.
If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE

Linearity
the method is linear up to 440 U/l.
If a ΔA/min of 0.200 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
the limit of detection is 0.463 U/l.

Interferences
no interference was observed by the presence of:
hemoglobin ≤ 250 mg/dl
bilirubin ≤ 45 mg/dl
lipids ≤ 500 mg/dl

Precision
intra-assay (n=10) mean (U/l) SD (U/l) CV%
sample 1 46.19 0.31 0.67
sample 2 137.25 0.92 0.67

inter-assay (n=20) mean (U/l) SD (U/l) CV%
sample 1 46.18 2.04 4.41
sample 2 137.76 6.30 4.57

Methods comparison
a comparison between Chema and a commercially available product gave the following results:

GOT Chema = x
GOT competitor = y
n = 83
y = 1.003x - 0.560 U/l r² = 0.990

WASTE DISPOSAL

This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.

GOT/AST FL IFCC	
GO 2H400	8 x 40 + 4 x 20 ml
GO 6U420	6 x 56 + 6 x 14 ml

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