# CRP CALIBRATOR

5 x 1 ml

REF RP CA05 CH

LOT -----

 -----

componente component composant componente	valore di calibrazione / calibration value / valeur de calibration / valor de calibración					unità unit unité unidad
	LEVEL 1	LEVEL 2	LEVEL 3	LEVEL 4	LEVEL 5	
Proteina C-reattiva C-Reactive Protein Protéine C-réactive Proteína C-reactiva	----	----	----	----	----	mg/l

## Italiano

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#### Uso previsto

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#### Bibliografia

"The certification of the mass concentration of C-Reactive Protein in human serum - Certified Reference Material ERM-DA474/IFCC", European Community report, 2011

## English

### HUMAN BASED CALIBRATOR SET

#### Intended use

CRP Calibrator is composed of a set of sera, aimed at the preparation of calibration curves for quantitative determination of C-Reactive Protein (CRP) in human serum, through immunoturbidimetric technique.

#### Principle of the method

C-Reactive Protein (CRP) selectively reacts with an anti-CRP antibody and forms an immunocomplex. The produced turbidity is proportional to the concentration of CRP in the sample, and can be measured at the wavelength of 340 nm.

#### Package components

Calibrators (liquid) 5 x 1 ml

#### Composition

CRP Calibrator is a liquid product. A dilution of defibrinated human plasma and recombinant human C-reactive protein with phosphate buffered saline. Containing stabilizers and sodium azide (< 0.1%).

#### Assigned values

CRP calibrator values are assigned to a reference material (ERM-DA474/IFCC) from the International Federation of Clinical Chemistry (IFCC).

#### Storage and stability

Store at 2-8°C. Do not freeze.  
After first opening the container, the serum can be used for 6 weeks if stored tightly closed 2-8°C after use.

#### Test procedure

Use supplied materials as provided.  
Additional materials required: automatic micropipette, distilled/deionized water, general laboratory equipment.

#### Precautions and warnings

For in vitro diagnostic use only.  
Exercise the normal precautions required for handling all laboratory reagents. Each individual donation intended for use in manufacture of protein standard serum was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required test. Only donations with negative findings were used for its manufacture. However, because no test method can offer complete assurance that products derived from human sources will not transmit infectious agents, it is recommended that this product is handled with the same precautions used for patient specimens.

#### Disposal

This product is made to be used in professional laboratories. Please consult local regulations for a correct waste disposal.  
P501: Dispose of contents according to national/international regulations.

#### References

"The certification of the mass concentration of C-Reactive Protein in human serum - Certified Reference Material ERM-DA474/IFCC", European Community report, 2011

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