# REFERENCE P MULTICALIBRATOR

5 x 1 ml

REF CP 0005 CH

LOT -----

-----

Français

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#### Bibliografia

"Certification of proteins in the human serum - Certified Reference Material ERM-DA470k/IFCC", European Community report, 2008

## English

### HUMAN BASED CALIBRATOR SET

#### Intended use

Reference P Multicalibrator is composed of a set of sera, aimed at the preparation of calibration curves for quantitative determination of proteins in human serum, through immunochemical techniques.

#### Package components

Calibrators (liquid) 5 x 1 ml

#### Composition

Reference P Multicalibrator is a liquid, defibrinated and delipidated human plasma, containing stabilizers, sodium azide (< 0.1%) as preservative and the following human proteins: albumin,  $\alpha$ 1-antitrypsin,  $\alpha$ 2-macroglobulin,  $\alpha$ 1-acid glycoprotein, ceruloplasmin, haptoglobin, transferrin, IgA, IgG, IgM, complement C3, complement C4, antithrombin III, C1 esterase inhibitor, kappa light chain, lambda light chain, prealbumin.

#### Assigned values

Protein values have been assigned to ERM-DA470k/IFCC, a Reference Preparation for proteins in Human Serum (RPPHS). RPPHS is a certified reference material developed by the College of American Pathologists (CAP), the Bureau Communautaire de Reference and the International Federation of Clinical Chemists (IFCC).

#### Storage and stability

Store at 2-8°C.

The calibrators are ready for use; after vial opening, sera are stable for 6 weeks if stored tightly closed at 2-8°C. Do not freeze.

#### Test procedure

Use supplied materials as provided.

Additional materials required: automatic micropipette, general laboratory equipment.

#### Precautions and warnings

For in vitro diagnostic use only.

Exercise the normal precautions required for handling all laboratory reagents. This product has been prepared exclusively from the blood of donors tested individually and shown by FDA approved methods to be free from HbsAg and antibodies to HCV and HIV. However, because no test method can offer complete assurance that products derived from human sources will not transmit infectious agents, it is recommended that this product be handled with the same precautions used for patient specimens.

#### Disposal

This product is made to be used in professional laboratories. Please consult local regulations for a correct waste disposal.

P501: Dispose of contents according to national/international regulations.

#### References

"Certification of proteins in the human serum - Certified Reference Material ERM-DA470k/IFCC", European Community report, 2008