**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ  
HITACHI 911/912**

TEST:	<b>CHOL</b>	
APP. CODE:	<b>263</b>	
WAVELENGTH (Sec/Pri):	<b>700 - 505</b>	
ASSAY:	<b>1-POINT</b>	<b>TIME: 10 POINT: 20</b>
SAMPLE VOL:	NORMAL: <b>3</b> DECREASE: <b>2</b> INCREASE: <b>5</b>	
	R1 VOLUME: <b>300</b> R2 VOLUME: <b>0</b> R3 VOLUME: <b>0</b> R4 VOLUME: <b>0</b>	
ABS LIMIT:	<b>32000 - INC</b>	
PROZONE LIMIT:	<b>0 - UPPER</b>	
CALIB METHOD:	<b>LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)</b>	
SD LIMIT:	<b>0.250</b>	
DUPLICATE LIMIT:	<b>3%</b>	
ST. 1 CONC:	<b>0.0</b>	
EXPECTED VALUE:	<b>100 - 200</b>	
UNIT:	<b>mg/dl</b>	
INSTR. FACTOR (y=ax+b):	a=1	b=0

**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ  
OLYMPUS AU 400/480/600/640/680/2700 (Test code 865)**

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**INTENDED USE**  
Reagent for quantitative in vitro determination of cholesterol in biological fluids.

**PRINCIPLE OF THE METHOD**  
All cholesterol esters present in plasma are hydrolyzed quantitatively into free cholesterol and fatty acids by cholesterol esterase. In the presence of oxygen, free cholesterol is then oxidized by cholesterol oxidase to colestene-4-ene-3-one and H<sub>2</sub>O<sub>2</sub>. The H<sub>2</sub>O<sub>2</sub> reacts with p-chlorophenol and 4-aminoantipyrine in the presence of peroxidase to form a quinoneimine dye. The intensity of color formed is proportional to the cholesterol concentration and can be measured photometrically between 480 and 520 nm.

**KIT COMPONENTS**  
**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources.

**CHOL R1 2H500 10 x 50 ml (liquid) white cap  
6U448 8 x 56 ml (liquid) white cap**

Composition: Good's buffer pH 7.20, sodium cholate 8 mM, CHE  $\geq$  400 U/l, CHOD  $\geq$  200 U/l, POD  $\geq$  500 U/l, 4-AAP 0.6 mM, 4-chlorophenol 2 mM.

Store all components at 2-8°C.

**REAGENT PREPARATION**  
Use separate reagent ready to use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.  
Caution: keep well refrigerated.

**CHOL R1: Warning.** Causes serious eye irritation (H319). Causes skin irritation (H315). Wear protective gloves. Eye protection (P280). IF ON SKIN: Wash with plenty of water (P302+P352). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). If eye irritation persists: get medical advice (P337+P313).

N-acetylcysteine (NAC), metamilozolo and acetaminophen may cause interference in the Trinder reaction.<sup>(1,2)</sup>  
To avoid interference, the blood withdrawal should be performed before drug administration.

**SPECIMEN**  
Serum, plasma EDTA.  
Sample is stable 3 days at 2-8°C and 1 month at -20°C.

<b>EXPECTED VALUES</b>	
desirable:	140 - 200 mg/dl
borderline/high risk:	200 - 240 mg/dl
high risk:	> 240 mg/dl

Each laboratory should establish appropriate reference intervals related to its population.

**QUALITY CONTROL AND CALIBRATION**

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:  
**QUANTINORM CHEMA**  
with normal or close to normal control values  
**QUANTIPATH CHEMA**  
with pathological control values.  
If required, a multiparametric, human based calibrator is available:  
**AUTOCAL H**

Please contact Customer Care for further information.

**TEST PERFORMANCE**

**Linearity**  
the method is linear up to 700 mg/dl.  
If the limit value is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**  
the limit of detection is 1 mg/dl.

**Interferences**  
no interference was observed by the presence of:  
hemoglobin  $\leq$  500 mg/dl  
bilirubin  $\leq$  15 mg/dl  
lipids  $\leq$  850 mg/dl

<b>Precision</b>				
intra-assay (n=10)	mean (mg/dl)	SD (mg/dl)	CV%	
sample 1	101.50	1.84	1.80	
sample 2	176.20	2.74	1.60	
inter-assay (n=20)	mean (mg/dl)	SD (mg/dl)	CV%	
sample 1	100.99	2.11	2.10	
sample 2	176.51	2.23	1.30	

**Methods comparison**  
a comparison between Chema and a commercially available product gave the following results:

$$\begin{aligned} \text{Cholesterol FL Chema} &= x \\ \text{Cholesterol competitor} &= y \\ n &= 100 \\ y &= 0.979x - 1.71 \text{ mg/dl} \quad r^2 = 0.995 \end{aligned}$$

**WASTE DISPOSAL**  
This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.