APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ <b>HITACHI 911/912</b>	
TEST:	<b>LIP</b>
APP. CODE:	<b>330</b>
WAVELENGTH (Sec/Pri):	<b>700 - 570</b>
ASSAY:	<b>RATE-A</b> <span style="float: right;"><i>TIME: 10</i> <i>POINT: 21 - 28</i> <i>DILUENT: water</i></span>
SAMPLE VOL:	NORMAL: <b>4</b> DECREASE: <b>3</b> INCREASE: <b>5</b>
	R1 VOLUME: <b>200</b> <span style="float: right;"><i>DILUENT: 5</i></span> R2 VOLUME: <b>0</b> R3 VOLUME: <b>50</b> <span style="float: right;"><i>DILUENT: 5</i></span> R4 VOLUME: <b>0</b>
ABS LIMIT:	<b>32000 - INC</b>
PROZONE LIMIT:	<b>0 - UPPER</b>
CALIB METHOD:	<b>LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)</b>
SD LIMIT:	<b>0.250</b>
DUPLICATE LIMIT:	<b>3%</b>
ST. 1 CONC:	<b>0.00</b>
EXPECTED VALUE:	<b>0 - 60</b>
UNIT:	<b>U/l</b>
INSTR. FACTOR (y=ax+b):	<b>a=1 b=0</b>

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ <b>OLYMPUS AU 400/480/600/640/680/2700 (Test code 879)</b>	
TEST NAME:	<b>LIP</b>
SAMPLE:	Volume <b>4 µl</b> <span style="float: right;">Dilution <b>0 µl</b></span>
REAGENTS:	R1 Volume <b>200 µl</b> <span style="float: right;">Dilution <b>0 µl</b></span> R2 Volume <b>50 µl</b> <span style="float: right;">Dilution <b>0 µl</b></span>
WAVELENGTH:	Pri. <b>570</b> Sec. <b>700</b>
METHOD:	<b>RATE</b>
REACTION SLOPE:	<b>+</b>
MEASURING POINT 1:	First <b>16</b> Last <b>24</b>
MEASURING POINT 2:	First Last
REAGENT OD LIMIT:	First L <b>0.01</b> First H <b>1.0</b> Last L <b>0.01</b> Last H <b>1.0</b>
DYNAMIC RANGE:	L <b>1</b> H <b>300</b>
CORRELATION FACTOR:	A <b>1</b> B <b>0</b>
UNIT:	<b>U/l</b>
CALIBRATION TYPE:	<b>AB</b>
FORMULA:	<b>Y = AX + B</b>

 Chema Diagnostica  
Via Campania 2/4  
60030 Monsano (AN) - ITALY - EU  
phone +39 0731 605064  
fax +39 0731 605672  
e-mail: mail@chema.com  
website: http://www.chema.com

ITALIANO rev. 12/07/2021

LIPASI FL	
LP 2H100	4 x 20 + 2 x 10 ml
LP 6U140	2 x 56 + 2 x 14 ml

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**INTENDED USE**  
Reagent for quantitative in vitro determination of lipase in biological fluids.

**PRINCIPLE OF THE METHOD**  
The colorimetric substrate 1,2-O-Dilauryl-rac-glycerol-3-glutaric acid-(6'-methyl-resorufin)-ester is cleaved by pancreatic lipase and the resulting dicarboxylic acid ester is hydrolysed under the alkaline test conditions to yield the chromophore methylresorufin. The kinetic of colour formation at 580 nm is monitored and it is proportional to lipase activity in sample.

**KIT COMPONENTS**  
**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label. Keep away from direct light sources.

**LIP R1** **2H100: 4 x 20 ml (liquid) white cap**  
**6U140: 2 x 56 ml (liquid) white cap**

Composition: Good's Buffer pH 8.0, colipase ≥ 1 mg/l, desoxycholate ≥ 1.0 mM, taurodesoxycholate ≥ 1.0 mM, calcium ions ≥ 1 mM, detergent and preservative.

**LIP R2** **2H100: 2 x 10 ml (liquid) red cap**  
**6U140: 2 x 14 ml (liquid) red cap**

Composition: Tartrate buffer pH 4.0, lipase substrate ≥ 0.1 mM, stabilizer and preservative.

Store all components at 2-8°C and do not freeze.

**REAGENT PREPARATION**  
Use separate reagents ready to use.  
Stability: up to expiration date on labels at 2-8°C;  
Stability since first opening of vials: use preferably within 60 days at 2-8°C.  
Caution: reagent R2 is a microemulsion. Therefore, a slight apparent precipitation could occur, showing a light red deposit on the bottom of vial. It is a normal behaviour and it is recommended to resuspend solution before analysis, with a mild shaking.

**PRECAUTIONS**  
**LIP R1:** It is not classified as hazardous.

**LIP R2: Danger.** Causes serious eye damage (H318).  
 Wear protective gloves. Eye protection (P280).  
IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). Immediately call a doctor (P310).

Some commercial reagents for triglycerides, HDL and LDL determination could contain microbial lipases, whose could stick on surface of instrument plastic cuvettes. It is recommended to program a "wash" procedure before lipase determination, if a contamination is suspected.  
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

**SPECIMEN**  
Serum, plasma heparinate.  
The lipase activity is stable for 7 days in samples stored at 2-8 ° C.

**EXPECTED VALUES**  
normal subjects: ≤ 60 U/l (methylresorufin 37°C)

Each laboratory should establish appropriate reference intervals related to its population.

**QUALITY CONTROL AND CALIBRATION**  
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:  
**QUANTINORM CHEMA - MULTINORM CHEMA** with normal or close to normal control values  
**QUANTIPATH CHEMA - MULTIPATH CHEMA** with pathological control values.  
If required, a multiparametric, human based calibrator is available:  
**AUTOCAL H**

Please contact Customer Care for further information.

**TEST PERFORMANCE**  
**Linearity**  
The method is linear up to 300 U/l.  
If the limit value is exceeded, it is suggested to dilute sample 1+1 with saline solution (9 g/l) and to repeat the test, multiplying the result by 2.

**Sensitivity/limit of detection (LOD)**  
The limit of detection is 1 U/l.

**Interferences**  
No interference was observed by the presence of:  
ascorbic acid ≤ 50 mg/dl  
hemoglobin ≤ 400 mg/dl  
bilirubin ≤ 50 mg/dl  
intralipid ≤ 1000 mg/dl

**Precision**  
intra-assay (n=10) mean (U/l) SD (U/l) CV%  
sample 1 49.9 0.65 1.30  
sample 2 110.5 1.69 1.53

inter-assay (n=20) mean (U/l) SD (U/l) CV%  
sample 1 50.0 1.43 2.87  
sample 2 110.9 3.91 3.53

**Methods comparison**  
A comparison between Chema and a commercially available product gave the following results:

Lipase Chema = y  
Lipase competitor = x  
n = 94  
 $y = 0.93x + 5.50$   $r^2 = 0.99$

**WASTE DISPOSAL**  
This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.