**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
HITACHI 911/912**

TEST:	ALB	
APP. CODE:	332	
WAVELENGTH (Sec/Pri):	700 - 600	
ASSAY:	1-POINT	<i>TIME: 10 POINT: 31</i>
SAMPLE VOL:	NORMAL: 2 DECREASE: 2 INCREASE: 5	
	R1 VOLUME: 300 R2 VOLUME: 0 R3 VOLUME: 0 R4 VOLUME: 0	
ABS LIMIT:	32000 - INC	
PROZONE LIMIT:	0 - UPPER	
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)	
SD LIMIT:	0.250	
DUPLICATE LIMIT:	3%	
ST. 1 CONC:	0.0	
EXPECTED VALUE:	0.0 - 5.5	
UNIT:	g/dl	
INSTR. FACTOR (y=ax+b):	a=1 b=0	

**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
OLYMPUS AU 400/480/600/640/680/2700 (Test code 857)**

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INTENDED USE
Reagent for quantitative in vitro determination of albumin in biological fluids.

PRINCIPLE OF THE METHOD
Albumin and BCG are followed to bind at pH 4.2, and absorption of the BCG-albumin complex is determined spectrophotometrically at 628 nm. At pH 4.2, albumin acts as a cation to bind the anionic dye.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

ALB R1 2H500: 10 x 50 ml (liquid) white cap
6U448: 8 x 56 ml (liquid) white cap

Composition: succinate buffer 100 mM pH 4.2, bromochresol green 0.2 mM, surfactant.
Store all components at 2-8°C.

REAGENT PREPARATION
Use reagent ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C.

PRECAUTIONS
Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SPECIMEN
Serum (preferred), plasma (heparinate or EDTA). Venostasis should be avoided in specimen collection because hemoconcentration increases the apparent concentrations of albumin and other plasma proteins.

EXPECTED VALUES

Men	4.2 - 5.5 g/dl
Women	3.7 - 5.3 g/dl

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:
QUANTINORM CHEMA with normal or close to normal control values
QUANTIPATH CHEMA with pathological control values.
If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE
Linearity
the method is linear up to 6 g/dl.
If the limit value is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
the limit of detection is 0.01 g/dl.

Interferences
no interference was observed by the presence of:
hemoglobin ≤ 350 mg/dl
bilirubin ≤ 27 mg/dl
lipids ≤ 850 mg/dl

Precision

intra-assay (n=10)	mean (g/dl)	SD (g/dl)	CV%
sample 1	3.37	0.04	1.10
sample 2	3.34	0.04	1.30

inter-assay (n=20)	mean (g/dl)	SD (g/dl)	CV%
sample 1	3.36	0.04	1.00
sample 2	3.35	0.07	2.00

Methods comparison
a comparison between Chema and a commercially available product gave the following results:

$$y = 1.009x - 0.195 \text{ g/dl} \quad r^2 = 0.956$$

WASTE DISPOSAL
This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.