APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	FEFZ
APP. CODE:	341
WAVELENGTH (Sec/Pri):	700 - 570
ASSAY:	2 POINT END <i>TIME: 10</i> <i>POINT: 16 - 31</i> <i>DILUENT: water</i>
SAMPLE VOL:	NORMAL: 50 DECREASE: 40 INCREASE: 60
	R1 VOLUME: 200 <i>DILUENT: 5</i> R2 VOLUME: 0 R3 VOLUME: 50 <i>DILUENT: 5</i> R4 VOLUME: 0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.00
EXPECTED VALUE:	30 - 160
UNIT:	µg/dl
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 871)	
TEST NAME:	FEFZ
SAMPLE:	Volume 40 µl Dilution 0 µl
REAGENTS:	R1 Volume 160 µl Dilution 0 µl R2 Volume 40 µl Dilution 0 µl
WAVELENGTH:	Pri. 570 Sec. 700
METHOD:	END
REACTION SLOPE:	+
MEASURING POINT 1:	First 0 Last 27
MEASURING POINT 2:	First 0 Last 10
REAGENT OD LIMIT:	First L -0.1 First H 0.5 Last L -0.1 Last H 0.5
DYNAMIC RANGE:	L 25 H 1000
CORRELATION FACTOR:	A 1 B 0
UNIT:	µg/dl
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

ITALIANO

rev. 11/06/2026

FERRO FZ	
FE 2H100	4 x 20 + 2 x 10 ml
FE 6U280	4 x 56 + 4 x 14 ml

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INTENDED USE
Reagent for quantitative in vitro determination of iron in biological fluids.

PRINCIPLE OF THE METHOD
Serum iron bound to transferrin is released in acidic environment. Fe(III) ions are then reduced to Fe(II), which reacts with ferrozine to give a violet colored complex. The absorbance measured at 560 nm is directly proportional to the amount of iron in the sample.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

FE FZ R1 2H100 4 x 20 ml (liquid) white cap
6U280 4 x 56 ml (liquid) white cap

Composition: acetate buffer 500 mM pH 4.50, thiourea ≥ 50 mM, guanidine hydrochloride ≥ 100 mM, surfactant.

FE FZ R2A 2H100 2 x 10 ml (liquid) red cap
6U280 4 x 14 ml (liquid) red cap
Composition: ferrozine 6 mM.

FE FZ R2B 2H100 2 vials powder for 10 ml
6U280 4 vials powder for 14 ml

Composition: sodium ascorbate ≥ 50 mM.

Store all components at 2-8°C.

REAGENT PREPARATION
Reagent R1: ready to use.
Reagent R2: add all the content of reagent R2B to reagent R2A and let to stay 20 minutes, mixing occasionally by inversion. Do not shake. Stable 90 days at 2-8°C.
Caution: keep well closed and refrigerated.
Stability of unmixed reagents:
up to expiration date on labels at 2-8°C.
Stability of unmixed reagents since first opening of vials: preferably within 60 days at 2-8°C.

PRECAUTIONS
FE FZ R1: Danger. Causes serious eye damage (H318). Causes skin irritation (H315). Harmful to aquatic life with long lasting effects (H412).
 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). Wear protective gloves / eye protection / face protection (P280). Immediately call a POISON CENTER / doctor (P310). Wash with water thoroughly after handling (P264). Avoid release to the environment (P273).

FE FZ R2A: It is not classified as hazardous.

FE FZ R2B: It is not classified as hazardous.

SPECIMEN
Serum, plasma heparinate.
Samples are stable 7 days at 15-25°C, 3 weeks at 2-8°C and several months at -20°C.
Separate serum/plasma from clot within 1 hour.
Anticoagulants as EDTA or oxalate could yield too low recovery values.

EXPECTED VALUES		
men	59 - 158 µg/dl	(10.6 - 28.3 µmol/l)
women	37 - 145 µg/dl	(6.60 - 26.0 µmol/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:
QUANTINORM CHEMA - MULTINORM CHEMA with normal or close to normal control values
QUANTIPATH CHEMA - MULTIPATH CHEMA with pathological control values.

Please contact Customer Care for further information.

TEST PERFORMANCE

Linearity
the method is linear up to 1000 µg/dl. If the limit value is exceeded, it is suggested to dilute sample 1+9 with distilled water and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
the limit of detection is 25 µg/dl.

Interferences
no interference was observed by the presence of:
hemoglobin interferes
bilirubin ≤ 19 mg/dl
lipids ≤ 1000 mg/dl

Precision	mean (µg/dl)	SD (µg/dl)	CV%
intra-assay (n=10)			
sample 1	106.41	2.12	1.99
sample 2	178.48	1.54	0.86
inter-assay (n=14)			
sample 1	107.69	6.65	6.20
sample 2	179.15	4.65	2.60

Methods comparison
a comparison between Chema and a commercially available product gave the following results:

$$y = 0.947x + 0.387 \mu\text{g/dl} \quad r^2 = 0.973$$

Iron FZ Chema = x
Iron competitor = y
n = 100

WASTE DISPOSAL

This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.

1/4