APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	GPT
APP. CODE:	328
WAVELENGTH (Sec/Pri):	700 - 340
ASSAY:	RATE-A
	<i>TIME: 10 POINT: 20-30 DILUENT: water</i>
SAMPLE VOL:	NORMAL: 25 DECREASE: 20 INCREASE: 30
	R1 VOLUME: 200 <i>DILUENT: 5</i> R2 VOLUME: 0 R3 VOLUME: 50 <i>DILUENT: 5</i> R4 VOLUME: 0
ABS LIMIT:	6500 - DEC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.0
EXPECTED VALUE:	0.0 - 45.0
UNIT:	U/l
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 854)	
TEST NAME:	GPT
SAMPLE:	Volume 25 µl Dilution 0 µl
REAGENTS:	R1 Volume 200 µl Dilution 0 µl R2 Volume 50 µl Dilution 0 µl
WAVELENGTH:	Pri. 340 Sec. 700
METHOD:	RATE
REACTION SLOPE:	-
MEASURING POINT 1:	First 14 Last 23
MEASURING POINT 2:	First Last
REAGENT OD LIMIT:	First L 0.8 First H 2.0 Last L 0.8 Last H 2.0
DYNAMIC RANGE:	L 0.2 H 440
CORRELATION FACTOR:	A 1 B 0
LINEARITY LIMIT:	15%
UNIT:	U/l
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

ITALIANO rev. 26/09/2016

GPT/ALT FL IFCC	
GP 2H400	8 x 40 + 4 x 20 ml
GP 6U420	6 x 56 + 6 x 14 ml

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INTENDED USE
Reagent for quantitative in vitro determination of GPT in biological fluids.

PRINCIPLE OF THE METHOD
The enzyme alanine aminotransferase (EC 2.6.1.2; L-Alanine:2-Oxoglutarate Aminotransferase, ALT or A1aAT; Glutamate Pyruvate Transaminase, GPT) catalyzes the transaminase reaction between L-Alanine and 2-Oxoglutarate. The pyruvate formed, is reduced to lactate in the presence of LDH. As the reactions proceed, NADH is oxidized to NAD. The disappearance of NADH per unit time is followed by measuring the decrease in absorbance at 340 nm. The present method has been made according to IFCC (2002).

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

GPT R1	2H400: 8 x 40 ml (liquid) white cap 6U420: 6 x 56 ml (liquid) white cap
GPT R2	2H400: 4 x 20 ml (liquid) red cap 6U420: 6 x 14 ml (liquid) red cap

Composition in the test: Tris buffer 100 mM pH 7.15, L-Alanine 500 mM, 2-Oxoglutarate 15 mM, NADH 0.18 mM, LDH ≥ 1700 U/l.

Store all components at 2-8°C.

REAGENT PREPARATION
Use separate reagent ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.
Caution: keep well refrigerated.

PRECAUTIONS
Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SPECIMEN
Serum (preferred). Plasma is not recommended.
Collect blood with a minimum of venous stasis.
GPT is stable up to 4 days at 2-8°C or 1 month at -20°C.

EXPECTED VALUES	
Men:	< 45 U/l (< 0.74 µkat/l)
Women:	< 34 U/l (< 0.56 µkat/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

QUANTINORM CHEMA
with normal or close to normal control values
QUANTIPATH CHEMA
with pathological control values.
If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE

Linearity
the method is linear up to 440 U/l.
If a ΔA/min of 0.200 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
the limit of detection is 0.169 U/l.

Interferences
no interference was observed by the presence of:
hemoglobin ≤ 500 mg/dl
bilirubin ≤ 45 mg/dl
lipids ≤ 450 mg/dl

PRECISION			
intra-assay (n=10)	mean (U/l)	SD (U/l)	CV%
sample 1	49.29	0.35	0.71
sample 2	132.15	0.57	0.43

INTER-SERIES (n=20)			
sample 1	mean (U/l)	SD (U/l)	CV%
sample 2	49.31	1.66	3.37
	132.85	4.28	3.22

Methods comparison
a comparison between Chema and a commercially available product gave the following results:

$$y = 0.992x - 0.299 \quad r^2 = 0.999$$

WASTE DISPOSAL
This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.