APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	LIP
APP. CODE:	330
WAVELENGTH (Sec/Pri):	700 - 570
ASSAY:	RATE-A <i>TIME: 10</i> <i>POINT: 21 - 28</i> <i>DILUENT: water</i>
SAMPLE VOL:	NORMAL: 4 DECREASE: 3 INCREASE: 5
	R1 VOLUME: 200 <i>DILUENT: 5</i> R2 VOLUME: 0 R3 VOLUME: 50 <i>DILUENT: 5</i> R4 VOLUME: 0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (<i>POINT: 2 - SPAN: 2 - WEIGHT: 0</i>)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.00
EXPECTED VALUE:	0 - 60
UNIT:	U/l
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 879)	
TEST NAME:	LIP
SAMPLE:	Volume 4 µl Dilution 0 µl
REAGENTS:	R1 Volume 200 µl Dilution 0 µl R2 Volume 50 µl Dilution 0 µl
WAVELENGTH:	Pri. 570 Sec. 700
METHOD:	RATE
REACTION SLOPE:	+
MEASURING POINT 1:	First 16 Last 24
MEASURING POINT 2:	First Last
REAGENT OD LIMIT:	First L 0.01 First H 1.0 Last L 0.01 Last H 1.0
DYNAMIC RANGE:	L 1 H 300
CORRELATION FACTOR:	A 1 B 0
UNIT:	U/l
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

ITALIANO rev. 12/04/2022

LIPASI FL	
LP 2H100	4 x 20 + 2 x 10 ml
LP 6U140	2 x 56 + 2 x 14 ml

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INTENDED USE
Reagent for quantitative in vitro determination of lipase in biological fluids.

PRINCIPLE OF THE METHOD
The colorimetric substrate 1,2-O-Dilauryl-rac-glycerol-3-glutaric acid-(6'-methyl-resorufin)-ester is cleaved by pancreatic lipase and the resulting dicarboxylic acid ester is hydrolysed under the alkaline test conditions to yield the chromophore methylresorufin. The kinetic of colour formation at 580 nm is monitored and it is proportional to lipase activity in sample.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label. Keep away from direct light sources.

LIP R1 **2H100: 4 x 20 ml (liquid) white cap**
6U140: 2 x 56 ml (liquid) white cap

Composition: Good's Buffer pH 8.0, colipase ≥ 1 mg/l, desoxycholate ≥ 1.0 mM, taurodesoxycholate ≥ 1.0 mM, calcium ions ≥ 1 mM, detergent and preservative.

LIP R2 **2H100: 2 x 10 ml (liquid) red cap**
6U140: 2 x 14 ml (liquid) red cap

Composition: Tartrate buffer pH 4.0, lipase substrate ≥ 0.1 mM, stabilizer and preservative.

Store all components at 2-8°C and do not freeze.

REAGENT PREPARATION
Use separate reagents ready to use.
Stability: up to expiration date on labels at 2-8°C;
Stability since first opening of vials: use preferably within 60 days at 2-8°C.
Caution: reagent R2 is a microemulsion. Therefore, a slight apparent precipitation could occur, showing a light red deposit on the bottom of vial. It is a normal behaviour and it is recommended to resuspend solution before analysis, with a mild shaking.

PRECAUTIONS
LIP R1: It is not classified as hazardous.

LIP R2: Danger. Causes serious eye damage (H318).
 Wear protective gloves. Eye protection (P280).
IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). Immediately call a doctor (P310).

Some commercial reagents for triglycerides, HDL and LDL determination could contain microbial lipases, whose could stick on surface of instrument plastic cuvettes. It is recommended to program a "wash" procedure before lipase determination, if a contamination is suspected.
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SPECIMEN
Serum, plasma heparinate.
The lipase activity is stable for 7 days in samples stored at 2-8 ° C.

EXPECTED VALUES
normal subjects: ≤ 60 U/l (methylresorufin 37°C)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:
QUANTINORM CHEMA - MULTINORM CHEMA with normal or close to normal control values
QUANTIPATH CHEMA - MULTIPATH CHEMA with pathological control values.
If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE
Linearity
The method is linear up to 300 U/l.
If the limit value is exceeded, it is suggested to dilute sample 1+1 with saline solution (9 g/l) and to repeat the test, multiplying the result by 2.

Sensitivity/limit of detection (LOD)
The limit of detection is 1 U/l.

Interferences
No interference was observed by the presence of:
ascorbic acid ≤ 50 mg/dl
hemoglobin ≤ 400 mg/dl
bilirubin ≤ 50 mg/dl
intralipid ≤ 1000 mg/dl

Precision
intra-assay (n=10) mean (U/l) SD (U/l) CV%
sample 1 49.9 0.65 1.30
sample 2 110.5 1.69 1.53

inter-assay (n=20) mean (U/l) SD (U/l) CV%
sample 1 50.0 1.43 2.87
sample 2 110.9 3.91 3.53

Methods comparison
A method comparison with the previous generation product gave the following results:

Chema current formulation = y
Chema previous formulation = x
n = 76
 $y = 1.017x - 1.452$ U/l $r^2 = 0.990$

WASTE DISPOSAL
This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.