APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912			
TEST: ALP DGKC			
APP. CODE: 342			
WAVELENGTH (Sec/Pri): 700 - 415			
ASSAY: RATE-A		TIME: 10 POINT: 20 - 30 DILUENT: water	
SAMPLE VOL:		NORMAL: 7 DECREASE: 5 INCREASE: 10	
		R1 VOLUME: 280 R2 VOLUME: 0 R3 VOLUME: 70 R4 VOLUME: 0	
		DILUENT: 5 DILUENT: 5	
ABS LIMIT: 32000 - INC			
PROZONE LIMIT: 0 - UPPER			
CALIB METHOD: LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)			
SD LIMIT: 0.250			
DUPLICATE LIMIT: 3%			
ST. 1 CONC: 0.0			
EXPECTED VALUE: 0 - 270			
UNIT: U/l			
INSTR. FACTOR (y=ax+b): a=1 b= 0			

APPLICAZIONE / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 872)			
TEST NAME: ALP DGKC			
SAMPLE: Volume 7 µl		Dilution 0 µl	
REAGENTS: R1 Volume 280 µl		Dilution 0 µl	
R2 Volume 70 µl		Dilution 0 µl	
WAVELENGTH: Pri. 410 Sec. 700			
METHOD: RATE			
REACTION SLOPE: +			
MEASURING POINT 1: First 14		Last 27	
MEASURING POINT 2: First		Last	
REAGENT OD LIMIT: First L -0.1		First H 1.0	
Last L -0.1		Last H 1.0	
DYNAMIC RANGE: L 7		H 3000	
CORRELATION FACTOR: A 1		B 0	
LINEARITY LIMIT: 15%			
UNIT: U/l			
CALIBRATION TYPE: AB			
FORMULA: Y = AX + B			

Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

FOSFATASI ALCALINA FL DGKC	
AL 2H100	4 x 20 + 2 x 10 ml
AL 2H400	8 x 40 + 4 x 20 ml
AL 6U420	6 x 56 + 6 x 14 ml

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INTENDED USE
Reagent for quantitative in vitro determination of alkaline phosphatase in biological fluids.
PRINCIPLE OF THE METHOD
The enzyme alkaline phosphatase (EC 3.1.3.1, ortho-phosphoric-monoester phosphohydrolase) hydrolizes the 4-NPP to release 4-nitrophenol, under alkaline conditions. The 4-nitrophenol formed is detected spectrophotometrically at 405 nm to give a measurement of alkaline phosphatase activity in the sample. The present method has been made according to DGKC.
KIT COMPONENTS
For in vitro diagnostic use only. The components of the kit are stable until expiration date on the label. Keep away from direct light sources.

ALP DGKC R1	2H100: 4 x 20 ml (liquid) white cap 2H400: 8 x 40 ml (liquid) white cap 6U420: 6 x 56 ml (liquid) white cap
ALP DGKC R2	2H100: 2 x 10 ml (liquid) red cap 2H400: 4 x 20 ml (liquid) red cap 6U420: 6 x 14 ml (liquid) red cap

Composition in the test: DEA buffer pH 9.8 1 M, MgCl₂ 0.5 mM, 4-Nitrophenilphosphate 10 mM.

Store all components at 2-8°C.
REAGENT PREPARATION
Use separate reagent ready to use. Stability: up to expiration date on labels at 2-8°C. Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-. Caution: keep well refrigerated.

PRECAUTIONS
ALP DGKC R1: Danger. May cause damage to organs through prolonged or repeated exposure (H373). Causes serious eye damage (H318). Causes skin irritation (H315). Wear protective gloves. Eye protection (P280). Do not breath vapours (P260). IF ON SKIN: Wash with plenty of water (P302+P352). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). Get medical advice if you feel unwell (P314).
 

ALP DGKC R2: It is not classified as hazardous.

SPECIMEN
Serum, plasma (heparinate only). Sera kept at room temperatures usually show a slight but real increase in activity, which varies from 1% over a 6-h period to 3 to 6% over a 1 to 4 days period. Even in sera stored at refrigerator temperature, activity increases slowly. In frozen sera, activity decreases but slowly recovers after thawing the serum. A similar enhancement of activity, but of greater magnitude, occurs with reconstituted lyophilized preparations, such as those available as control sera or calibrators. In reconstituted material the increases with storage at 4 and 20°C are about 10 and 30%, respectively. Enhancement of activity continues for several days, but at a decreasing rate. The cause of this phenomenon is not known but may be attributed to renaturation of partially denatured enzyme or to dissociation, on warming, of a phosphate-lipoprotein complex or a multimer of the enzyme that was formed in the freeze-drying process.
EXPECTED VALUES
Men: < 270 U/l (< 4.50 µkat/l) Women: < 240 U/l (< 4.00 µkat/l)
QUALITY CONTROL AND CALIBRATION
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available: QUANTINORM CHEMA with normal or close to normal control values QUANTIPATH CHEMA with pathological control values. If required, a multiparametric, human based calibrator is available: AUTOCAL H

EXPECTED VALUES
Men: < 270 U/l (< 4.50 µkat/l) Women: < 240 U/l (< 4.00 µkat/l)

QUALITY CONTROL AND CALIBRATION
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available: QUANTINORM CHEMA with normal or close to normal control values QUANTIPATH CHEMA with pathological control values. If required, a multiparametric, human based calibrator is available: AUTOCAL H

Please contact Customer Care for further informations.

TEST PERFORMANCE
Linearity the method is linear up to 3000 U/l. If a ΔA/min of 0.500 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
the limit of detection is 7 U/l.

Interferences no interference was observed by the presence of: hemoglobin bilirubin lipids	≤ 400 mg/dl ≤ 27 mg/dl ≤ 1000 mg/dl
---	---

Precision intra-assay (n=10) sample 1 sample 2	mean (U/l) 175.70 426.70	SD (U/l) 0.95 2.41	CV% 0.50 0.60
inter-assay (n=20) sample 1 sample 2	mean (U/l) 167.26 408.28	SD (U/l) 3.99 8.61	CV% 2.40 2.10

Methods comparison
a comparison between Chema and a commercially available product gave the following results:

ALP Chema = x ALP competitor = y n = 112	y = 0.96x - 2.17 U/l r² = 0.999
--	------------------------------------

WASTE DISPOSAL
This product is made to be used in professional laboratories. P501: Dispose of contents according to national/international regulations.

PHOSPHATASE ALCALINE FL DGKC		
AL 2H100	4 x 20 + 2 x 10 ml	
AL 2H400	8 x 40 + 4 x 20 ml	
AL 6U420	6 x 56 + 6 x 14 ml	

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