APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	AOX
APP. CODE:	331
WAVELENGTH (Sec/Pri):	700 - 546
ASSAY:	2 POINT END <i>TIME: 10</i> <i>POINT: 16 - 31</i> <i>DILUENT: water</i>
SAMPLE VOL:	NORMAL: 12 DECREASE: 10 INCREASE: 15
	R1 VOLUME: 240 <i>DILUENT: 5</i> R2 VOLUME: 0 R3 VOLUME: 60 <i>DILUENT: 5</i> R4 VOLUME: 0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.0
EXPECTED VALUE:	0.0 - 7.2
UNIT:	mg/dl
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 855)	
TEST NAME:	AOX
SAMPLE:	Volume 10 µl Dilution 0 µl
REAGENTS:	R1 Volume 200 µl Dilution 0 µl R2 Volume 50 µl Dilution 0 µl
WAVELENGTH:	Pri. 540 Sec. 700
METHOD:	END
REACTION SLOPE:	+
MEASURING POINT 1:	First 0 Last 27
MEASURING POINT 2:	First 0 Last 10
REAGENT OD LIMIT:	First L -0.1 First H 0.5 Last L -0.1 Last H 0.5
DYNAMIC RANGE:	L 0.05 H 35
CORRELATION FACTOR:	A 1 B 0
UNIT:	mg/dl
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

ITALIANO rev. 09/05/2019

ACIDO URICO AOX FL	
AX 2H200	4 x 40 + 2 x 20 ml
AX 6U420	6 x 56 + 6 x 14 ml

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INTENDED USE
Reagent for quantitative in vitro determination of uric acid in biological fluids.

PRINCIPLE OF THE METHOD
Uric acid in sample is oxidized to allantoin in presence of the enzyme uricase and H₂O₂ is generated. The H₂O₂ reacts with TOOS and 4-aminoantipyrine in the presence of peroxidase to form a violet dye. The intensity of color formed is proportional to the uric acid concentration and can be measured photometrically between 510 and 560 nm. Product contains ascorbate oxidase to eliminate ascorbic acid interference.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

UA AOX R1 2H200: 4 x 40 ml (liquid) white cap
6U420: 6 x 56 ml (liquid) white cap

Composizione: phosphate buffer pH 7.0 100 mM, TOOS 0.38 mM, ascorbate oxidase ≥ 1000 U/l, surfactant.

UA AOX R2 2H200: 2 x 20 ml (liquid) red cap
6U420: 6 x 14 ml (liquid) red cap

Composizione: Good buffer pH 7.7 50 mM, 4-aminoantipyrine 1.5 mM, uricase ≥ 450 U/l, POD ≥ 1000 u/l, surfactant.

Store all components at 2-8°C.

REAGENT PREPARATION
Use separate reagents ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C

PRECAUTIONS

UA AOX R1: Warning. Causes serious eye irritation (H319). Causes skin irritation (H315). Wear protective gloves. Eye protection (P280). IF ON SKIN: Wash with plenty of water (P302+P352). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). If eye irritation persists: get medical advice (P337+P313).

UA AOX R2: Warning. Causes serious eye irritation (H319). Causes skin irritation (H315). Wear protective gloves. Eye protection (P280). IF ON SKIN: Wash with plenty of water (P302+P352). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). If eye irritation persists: get medical advice (P337+P313).

N-acetylcysteine (NAC), metamizole and acetaminophen may cause interference in the Trinder reaction.^(1,2) To avoid interference, the blood withdrawal should be performed before drug administration.

SPECIMEN
Serum, plasma heparinate. Oxalate, citrate and fluoride could yield a small decrease of uric acid. Urine. Uric acid is stable 5 days at 4-25°C. Dilute urine sample 1:10 with deionized water.

EXPECTED VALUES
Serum/plasma samples:
Men: 3.5 - 7.2 mg/dl (0.21 - 0.42 mmol/l)
Women: 2.6 - 6.0 mg/dl (0.15 - 0.35 mmol/l)

24h urine:
250 - 750 mg/24h (1.50 - 4.50 mmol/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:
QUANTINORM CHEMA
with normal or close to normal control values
QUANTIPATH CHEMA
with pathological control values.
If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE
Linearity
the method is linear up to 35 mg/dl.
If the value is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
the limit of detection is 0.06 mg/dl.

Interferences
no interference was observed by the presence of:
hemoglobin ≤ 1000 mg/dl
bilirubin ≤ 29 mg/dl
lipids ≤ 970 mg/dl
ascorbic acid ≤ 50 mg/dl

Precision
intra-assay (n=10) mean (mg/dl) SD (mg/dl) CV%
sample 1 4.49 0.02 0.47
sample 2 12.04 0.06 0.49
inter-assay (n=21) mean (mg/dl) SD (mg/dl) CV%
sample 1 4.53 0.08 1.67
sample 2 12.01 0.24 2.00

Methods comparison
a comparison between Chema and a commercially available product gave the following results:

Uric acid AOX FL Chema = x
Uric acid competitor = y
n = 120

y = 0.882 x + 0.037 mg/dl r² = 0.99

WASTE DISPOSAL
This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.