**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
HITACHI 911/912**

TEST:	TRIG	
APP. CODE:	351	
WAVELENGTH (Sec/Pri):	700 - 546	
ASSAY:	1-POINT	TIME: 10 POINT: 20
SAMPLE VOL:	NORMAL: 3 DECREASE: 2 INCREASE: 5	
	R1 VOLUME: 300 R2 VOLUME: 0 R3 VOLUME: 0 R4 VOLUME: 0	
ABS LIMIT:	32000 - INC	
PROZONE LIMIT:	0 - UPPER	
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)	
SD LIMIT:	0.250	
DUPLICATE LIMIT:	3%	
ST. 1 CONC:	0.0	
EXPECTED VALUE:	0.0 - 200	
UNIT:	mg/dl	
INSTR. FACTOR (y=ax+b):	a=1	b=0

**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
OLYMPUS AU 400/480/600/640/680/2700 (Test code 883)**

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INTENDED USE
Reagent for quantitative in vitro determination of triglycerides in biological fluids.

PRINCIPLE OF THE METHOD
Triglycerides are hydrolyzed by lipoproteinlipase to produce glycerol and free fatty acids. The glycerol participates in a series of coupled enzymatic reactions, in which glycerol kinase / glycerol phosphate oxidase are involved and H₂O₂ is generated. H₂O₂ reacts with TOPS and 4-aminoantipyrine in the presence of peroxidase to form a quinoneimine dye. The intensity of color formed is proportional to the triglycerides concentration and can be measured photometrically at 546 nm.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label at 2-8°C.
Keep away from direct light sources.

**TRIG R1 2H500: 10 x 50 ml (liquid) white cap
6U448: 8 x 56 ml (liquid) white cap**

Composition: Good's buffer pH 6.80, ATP 2 mM, GK > 300 U/l, POD > 1000 U/l, LPL > 1000 U/l, GPO > 2000 U/l, TOPS 3 mM, 4-AAP 0.3 mM, , surfactants and stabilizers.

Store all components at 2-8°C.

REAGENT PREPARATION
Use separate reagent ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.
Caution: keep well refrigerated.

TRIG R1: Warning. Causes serious eye irritation (H319). Causes skin irritation (H315). Wear protective gloves. Eye protection (P280). IF ON SKIN: Wash with plenty of water (P302+P352). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). If eye irritation persists: get medical advice (P337+P313).

N-acetylcysteine (NAC), metamizole and acetaminophen may cause interference in the Trinder reaction.^(1,2)
To avoid interference, the blood withdrawal should be performed before drug administration.

SPECIMEN
Specimens should not be obtained for triglyceride determination unless the patient has been fasting for 10 to 14 h. Either serum or EDTA plasma can be used to determine triglycerides. When EDTA plasma is used, the plasma value is converted to the equivalent serum value by multiplying the plasma value by 1.03. Store specimens at 4°C before analysis. Specimens are stable at 4°C for 3 days, frozen at -20°C for two weeks, or frozen at -70°C for longer periods. Lipemic specimens may require warming to 37°C and vigorous mixing before analysis, especially if they have been frozen.

EXPECTED VALUES
desirable: < 200 mg/dl (2.26 mmol/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:
QUANTINORM CHEMA
with normal or close to normal control values
QUANTIPATH CHEMA
with pathological control values.
If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE
Linearity
the method is linear up to 1000 mg/dl.
If the value is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
the limit of detection is 0.69 mg/dl.

Interferences
no interference was observed by the presence of:
hemoglobin ≤ 150 mg/dl
bilirubin ≤ 18 mg/dl

Precision	intra-assay (n=10)	mean (mg/dl)	SD (mg/dl)	CV%
	sample 1	109.61	1.02	0.93
	sample 2	214.62	1.10	0.51

	inter-assay (n=20)	mean (mg/dl)	SD (mg/dl)	CV%
	sample 1	108.64	3.31	3.05
	sample 2	210.25	6.54	3.11

Methods comparison
a comparison between Chema TRIGLYCERIDES FL and a commercially available product gave the following results:

Triglycerides Chema = x	Triglycerides competitor = y	
n = 96		
y = 0.9993 x - 0.614 mg/dl	r ² = 0.995	

WASTE DISPOSAL
This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.