APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN <b>HITACHI 911/912</b>	
TEST:	<b>ASO</b>
APP. CODE:	<b>357</b>
WAVELENGTH (Sec/Pri):	<b>600</b>
ASSAY:	<b>2 POINT END</b> <i>TIME: 10</i> <i>POINT: 16 - 31</i> <i>DILUENT: water</i>
SAMPLE VOL:	NORMAL: <b>3</b> DECREASE: <b>2</b> INCREASE: <b>5</b>
	R1 VOLUME: <b>250</b> R2 VOLUME: <b>0</b> R3 VOLUME: <b>40</b> <i>DILUENT: 3</i> R4 VOLUME: <b>0</b>
ABS LIMIT:	<b>32000 - INC</b>
PROZONE LIMIT:	<b>0 - UPPER</b>
CALIB METHOD:	<b>LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)</b>
SD LIMIT:	<b>0.50</b>
DUPLICATE LIMIT:	<b>5%</b>
ST. 1 CONC:	<b>0.0</b>
EXPECTED VALUE:	<b>0.0 - 200</b>
UNIT:	<b>IU/ml</b>
INSTR. FACTOR (y=ax+b):	a=1    b= 0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN <b>OLYMPUS AU 400/480/600/640/680/2700 (Test code 892)</b>	
TEST NAME:	<b>ASO</b>
SAMPLE:	Volume <b>3</b> µl      Dilution <b>0</b> µl
REAGENTS:	R1 Volume <b>250</b> µl      Dilution <b>0</b> µl R2 Volume <b>40</b> µl      Dilution <b>0</b> µl
WAVELENGTH:	Pri. <b>600</b> Sec.
METHOD:	<b>END</b>
REACTION SLOPE:	<b>+</b>
MEASURING POINT 1:	First <b>0</b> Last <b>27</b>
MEASURING POINT 2:	First <b>0</b> Last <b>10</b>
REAGENT OD LIMIT:	First L <b>-0.1</b> First H <b>0.5</b> Last L <b>-0.1</b> Last H <b>0.5</b>
DYNAMIC RANGE:	L <b>10</b> H <b>400</b>
CORRELATION FACTOR:	A <b>1</b> B <b>0</b>
UNIT:	<b>IU/ml</b>
CALIBRATION TYPE:	<b>AB</b>
FORMULA:	<b>y = Ax + B</b>

 Chema Diagnostica  
Via Campania 2/4  
60030 Monsano (AN) - ITALY - EU  
phone +39 0731 605064  
fax +39 0731 605672  
e-mail: mail@chema.com  
website: http://www.chema.com

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**INTENDED USE**  
Reagent for quantitative in vitro determination of antistreptolysin O in biological fluids.

**PRINCIPLE OF THE METHOD**  
Antibodies of the sample selectively react with Streptolysin O coated to latex, thus producing particle agglutination. The produced turbidity is proportional to the concentration of ASO in the sample, and can be measured at the wavelength of 600 nm.

**KIT COMPONENTS**  
**For in vitro diagnostic use only.**  
The components of the kit are stable until expiration date on the label.  
Keep away from direct light sources.

**ASO R1 2H100 2 x 42 ml (liquid) white cap**  
**6U140 2 x 54 ml (liquid) white cap**

Composizione: Buffer pH 7.4, stabilizers and preservatives.

**ASO R2 2H100 1 x 14 ml (liquid) red cap**  
**6U140 2 x 9 ml (liquid) red cap**

Composizione: suspension of latex particles coated with Streptolysin O, stabilizers and preservatives.

Store all components at 2-8°C.

**REAGENT PREPARATION**  
Use separate reagents ready to use.  
Kindly shake R2 vial before use.  
Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: use preferably within 60 days at 2-8°C.

**PRECAUTIONS**  
Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.  
Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

**SPECIMEN**  
Serum, plasma.  
Keep specimens away from direct light sources.  
Samples are stable 2 days when stored at 2-8°C and 1 month at -20°C.

**EXPECTED VALUES**  
Adults < 200 IU/ml

Each laboratory should establish appropriate reference intervals related to its population.

**QUALITY CONTROL AND CALIBRATION**

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

**QUANTINORM CHEMA - MULTINORM CHEMA**  
with normal or close to normal control values,  
**QUANTIPATH CHEMA - MULTIPATH CHEMA**  
with pathological or close to pathological control values.

Please contact Customer Care for further information.

**TEST PERFORMANCE**  
**Linearity**  
the method is linear up to 400 IU/ml.  
If the value is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

**Hook effect**  
No Hook effect is observed with concentrations lower than 4100 IU/ml.

**Sensitivity/limit of detection**  
The limit of detection is 10 IU/ml.

**Interferences**  
No interference was observed by the presence of:  
hemoglobin ≤ 1000 mg/dl  
bilirubin ≤ 60 mg/dl  
lipids ≤ 2500 mg/dl  
rheumatoid factor ≤ 250 IU/ml

**Precision**  
intra-assay (n=10) mean (IU/ml) SD (IU/ml) CV%  
sample 1 74.7 1.56 2.09  
sample 2 173.1 2.84 1.64

inter-assay (n=20) mean (IU/ml) SD (IU/ml) CV%  
sample 1 75.7 3.44 4.54  
sample 2 172.7 4.58 2.65

**Methods comparison**  
A comparison between CHEMA and a commercially available product gave the following results:

ASO competitor = x  
ASO FL CHEMA = y  
n = 43  
 $y = 0.93x + 2.88 \text{ IU/ml}$        $r^2 = 0.97$

**WASTE DISPOSAL**  
This product is made to be used in professional laboratories.  
P501: Dispose of contents according to national/international regulations.

ASO FL	
AS 2H100	2 x 42 + 1 x 14 ml
AS 6U140	2 x 54 + 2 x 9 ml

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