**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
HITACHI 911/912**

TEST:	AMY	
APP. CODE:	334	
WAVELENGTH (Sec/Pri):	700 - 415	
ASSAY:	RATE-A	TIME: 10 POINT: 7 - 14
SAMPLE VOL:	NORMAL: 7 DECREASE: 5 INCREASE: 10	
	R1 VOLUME: 280 R2 VOLUME: 0 R3 VOLUME: 0 R4 VOLUME: 0	
ABS LIMIT:	32000 - INC	
PROZONE LIMIT:	0 - UPPER	
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)	
SD LIMIT:	0.250	
DUPLICATE LIMIT:	3%	
ST. 1 CONC:	0.00	
EXPECTED VALUE:	0 - 96	
UNIT:	U/l	
INSTR. FACTOR (y=ax+b):	a=1 b=0	

**APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ
OLYMPUS AU 400/480/600/640/680/2700 (Test code 859)**

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INTENDED USE

Reagent for quantitative in vitro determination of amylase in biological fluids.

PRINCIPLE OF THE METHOD

The enzyme α-amylase hydrolyzes the 2-chloro-4-nitrophenyl-α-D-maltotrioside (CNP-G3) to release 2-chloro-4-nitrophenol and form 2-chloro-4-nitrophenyl-α-D-maltoside (CNP-G2), maltotriose (G3) and glucose (G). The rate of formation of the 2-chloro-4-nitrophenol can be detected spectrophotometrically at 405 nm to give a measurement of α-amylase activity in the sample.

KIT COMPONENTS

For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

AMY R1 2H120: 6 x 20 ml (liquid) white cap
6U168: 3 x 56 ml (liquid) white cap

Composition: CNP-G3 2.3 mM, NaCl 350 mM, calcium acetate 6 mM, potassium thiocyanate 600 mM, Good's buffer pH 6.0 100 mM, stabilizers and non-reactive components.

Store all components at 2-8°C.

REAGENT PREPARATION

Use reagent ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C.

PRECAUTIONS

AMY R1: Danger. Causes serious eye damage (H318).
 Wear protective gloves. Eye protection (P280). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). Immediately call a doctor (P310).

SPECIMEN

Serum, plasma (heparinate only). Urine.
Amylase is stable in serum and plasma sample up to 2 months at 2-8°C.

EXPECTED VALUES

Serum - plasma: < 96 U/l (< 1.60 µkat/l)
Urine: < 480 U/l (< 8.00 µkat/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

QUANTINORM CHEMA
with normal or close to normal control values
QUANTIPATH CHEMA
with pathological control values.
If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE

Linearity
the method is linear up to 2000 U/l.
If a ΔA/min of 0.500 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
the limit of detection is 0.91 U/l.

Interferences
no interference was observed by the presence of:
hemoglobin ≤ 500 mg/dl
bilirubin ≤ 50 mg/dl
lipids ≤ 1200 mg/dl

Precision

intra-assay (n=10)	mean (U/l)	SD (U/l)	CV%
sample 1	67.89	0.97	1.42
sample 2	171.67	2.61	1.52

inter-assay (n=21)	mean (U/l)	SD (U/l)	CV%
sample 1	67.81	1.93	2.85
sample 2	175.16	4.92	2.81

Methods comparison
a comparison between Chema and a commercially available product gave the following results:

$$\begin{aligned} \text{Amylase Chema} &= x \\ \text{Amylase competitor} &= y \\ n &= 155 \\ y &= 1.071x - 0.54 \text{ U/l} \quad r^2 = 0.997 \end{aligned}$$

WASTE DISPOSAL

This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.