APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	EPS
APP. CODE:	333
WAVELENGTH (Sec/Pri):	700 - 415
ASSAY:	RATE-A <i>TIME: 10 POINT: 14 - 25 DILUENT: water</i>
SAMPLE VOL:	NORMAL: 10 DECREASE: 8 INCREASE: 12
R1 VOLUME:	240 <i>DILUENT: 5</i>
R2 VOLUME:	60 <i>DILUENT: 5</i>
R3 VOLUME:	0
R4 VOLUME:	0
ABS LIMIT:	32000 - INC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.00
EXPECTED VALUE:	0 - 100
UNIT:	U/l
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 858)	
TEST NAME:	EPS
SAMPLE:	Volume 10 µl Dilution 0 µl
REAGENTS:	R1 Volume 240 µl Dilution 0 µl R2 Volume 60 µl Dilution 0 µl
WAVELENGTH:	Pri. 410 Sec. 700
METHOD:	RATE
REACTION SLOPE:	+
MEASURING POINT 1:	First 17 Last 27
MEASURING POINT 2:	First Last
REAGENT OD LIMIT:	First L -0.1 First H 0.5 Last L -0.1 Last H 0.5
DYNAMIC RANGE:	L 6 H 1500
CORRELATION FACTOR:	A 1 B 0
LINEARITY LIMIT:	15%
UNIT:	U/l
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

 Chema Diagnostica
 Via Campania 2/4
 60030 Monsano (AN) - ITALY - EU
 phone +39 0731 605064
 fax +39 0731 605672
 e-mail: mail@chema.com
 website: http://www.chema.com

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INTENDED USE
 Reagent for quantitative in vitro determination of amylase in biological fluids.

PRINCIPLE OF THE METHOD
 The enzyme α-amylase (EC 3.2.1.1, 1,4 α-D-glucose glucanohydrolase) hydrolyzes the EPS to release several different fragments. The fragments so formed are completely hydrolyzed to 4-nitrophenol and glucose by α-glucosidase. The 4-nitrophenol formed is detected spectrophotometrically at 405 nm to give a measurement of α-amylase activity in the sample.
 The present method has been made according to IFCC.

KIT COMPONENTS
For in vitro diagnostic use only.
 The components of the kit are stable until expiration date on the label.
 Keep away from direct light sources.

**AMY EPS R1 2H100: 4 x 20 ml (liquid) white cap
 6U280: 4 x 56 ml (liquid) white cap**

**AMY EPS R2 EA 2H100: 2 x 10 ml (liquid) red cap
 EA 6U280: 4 x 14 ml (liquid) red cap**

Composition in the test: Hepes buffer 50 mM pH 7.10, NaCl 70 mM, calcium acetate 1.0 mM, EPS-G7 5.0 mM, α-glucosidase 6 kU/l.

Store all components at 2-8°C.

REAGENT PREPARATION
 Use reagent ready to use.
 Stability: up to expiration date on labels at 2-8°C.
 Stability since first opening of vials: preferably within 60 days at 2-8°C.

PRECAUTIONS
 Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.
 Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

SPECIMEN
 Serum, plasma (heparinate only), Urine.
 Amylase is stable in serum and plasma sample up to 2 months at 2-8°C.

EXPECTED VALUES

Serum - plasma:	28 - 100 U/l (0.47 - 1.67 µkat/l)
Random urine:	≤ 460 U/l (≤ 7.68 µkat/l)

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION

It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

QUANTINORM CHEMA
 with normal or close to normal control values
QUANTIPATH CHEMA
 with pathological control values.

If required, a multiparametric, human based calibrator is available:
AUTOCAL H

Please contact Customer Care for further information.

TEST PERFORMANCE

Linearity
 the method is linear up to 1500 U/l.
 If a ΔA/min of 0.500 is exceeded, it is suggested to dilute sample 1+9 with saline and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
 the limit of detection is 6 U/l.

Interferences
 no interference was observed by the presence of:
 hemoglobin ≤ 200 mg/dl
 bilirubin ≤ 48 mg/dl
 lipids interfer in low values

Precision

intra-assay (n=10)	mean (U/l)	SD (U/l)	CV%
sample 1	77.90	0.74	0.90
sample 2	194.80	1.99	1.00

inter-assay (n=20)

mean (U/l)	SD (U/l)	CV%
sample 1	75.77	1.90
sample 2	194.15	2.39

Methods comparison
 a comparison between Chema Amylase FL and a commercially available product gave the following results:

Amylase Chema = x	
Amylase competitor = y	
n = 108	
y = 1.067x + 5.21 U/l	r ² = 0.99

WASTE DISPOSAL
 This product is made to be used in professional laboratories.
 P501: Dispose of contents according to national/international regulations.