APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	<b>CKNAC</b>
APP. CODE:	<b>336</b>
WAVELENGTH (Sec/Pri):	<b>700 - 340</b>
ASSAY:	<b>RATE-A</b> <i>TIME: 10 POINT: 24 - 31 DILUENT: water</i>
SAMPLE VOL:	NORMAL: <b>10</b> DECREASE: <b>7</b> INCREASE: <b>12</b>
R1 VOLUME:	<b>200</b> <i>DILUENT: 5</i>
R2 VOLUME:	<b>0</b>
R3 VOLUME:	<b>50</b> <i>DILUENT: 5</i>
R4 VOLUME:	<b>0</b>
ABS LIMIT:	<b>32000 - INC</b>
PROZONE LIMIT:	<b>0 - UPPER</b>
CALIB METHOD:	<b>LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)</b>
SD LIMIT:	<b>0.250</b>
DUPLICATE LIMIT:	<b>3%</b>
ST. 1 CONC:	<b>0.00</b>
EXPECTED VALUE:	<b>0 - 200</b>
UNIT:	<b>U/I</b>
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 864)	
TEST NAME:	<b>CKNAC</b>
SAMPLE:	Volume <b>10 µl</b> Dilution <b>0 µl</b>
REAGENTS:	R1 Volume <b>200 µl</b> Dilution <b>0 µl</b> R2 Volume <b>50 µl</b> Dilution <b>0 µl</b>
WAVELENGTH:	Pri. <b>340</b> Sec. <b>700</b>
METHOD:	<b>RATE</b>
REACTION SLOPE:	<b>+</b>
MEASURING POINT 1:	First <b>18</b> Last <b>26</b>
MEASURING POINT 2:	First Last
REAGENT OD LIMIT:	First L <b>-0.1</b> First H <b>1.0</b> Last L <b>-0.1</b> Last H <b>1.0</b>
DYNAMIC RANGE:	L <b>1</b> H <b>2000</b>
CORRELATION FACTOR:	A <b>1</b> B <b>0</b>
LINEARITY LIMIT:	<b>15%</b>
UNIT:	<b>U/I</b>
CALIBRATION TYPE:	<b>AB</b>
FORMULA:	<b>Y = AX + B</b>

 Chema Diagnostica  
 Via Campania 2/4  
 60030 Monsano (AN) - ITALY - EU  
 phone +39 0731 605064  
 fax +39 0731 605672  
 e-mail: mail@chema.com  
 website: http://www.chema.com

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**INTENDED USE**  
 Reagent for quantitative in vitro determination of creatine kinase in biological fluids.

**PRINCIPLE OF THE METHOD**  
 Creatine kinase (EC 2.7.3.2; adenosine triphosphate: creatine N-phosphotransferase; CK) catalyzes the conversion of creatine phosphate and ADP to creatine and ATP. ATP and glucose are converted to ADP and glucose-6-phosphate by hexokinase. Glucose-6-phosphate dehydrogenase oxidizes glucose-6-phosphate to 6-phosphogluconate, reducing NADP to NADPH. The rate of conversion of NADP/NADPH, monitored at 340 nm, is proportional to CK activity. N-acetyl cysteine (NAC) is added as an activator of CK.

**KIT COMPONENTS**  
**For in vitro diagnostic use only.**  
 The components of the kit are stable until expiration date on the label.  
 Keep away from direct light sources.

**CK-NAC R1** 2H100: 4 x 20 ml (liquid) white cap  
 6U280: 4 x 56 ml (liquid) white cap

**CK-NAC R2** 2H100: 2 x 10 ml (liquid) red cap  
 6U280: 4 x 14 ml (liquid) red cap

Composition in the test: imidazole buffer 29 mM pH 6.50, creatine phosphate 30 mM, glucose 20 mM, N-acetyl-L-cysteine 20 mM, magnesium acetate 10 mM, EDTA 2 mM, ADP 2 mM, NADP 2 mM, AMP 5 mM, D(adenosine-5') pentaphosphate 12 µM, glucose-6-phosphate-dehydrogenase ≥ 3 kU/l, hexokinase ≥ 3 kU/l.

Store all components at 2-8°C.

**REAGENT PREPARATION**  
 Use separate reagent ready to use.  
 Stability: up to expiration date on labels at 2-8°C.  
 Stability since first opening of vials: preferably within 60 days at 2-8°C - away from light sources.  
 Caution: keep well refrigerated.

**PRECAUTIONS**  
 Reagent may contain some non-reactive and preservative components. It is suggested to handle carefully it, avoiding contact with skin and swallow.  
 Perform the test according to the general "Good Laboratory Practice" (GLP) guidelines.

**SPECIMEN**  
 Serum is the preferred specimen. Plasma containing heparin, EDTA, citrate, or fluoride may produce unpredictable reaction rates. CK activity in serum is unstable and is rapidly lost during storage. CK is inactivated both by bright daylight and by increasing specimen pH owing to loss of carbon dioxide; accordingly, specimens should be stored in the dark in tightly closed tubes. CK is susceptible to thermal denaturation; the degree of inactivation corresponds to the degree of temperature increase. Therefore, the serum specimen should be chilled to 4°C as rapidly as possible after collection. A slight degree of hemolysis can be tolerated because erythrocytes contain no CK activity. However, moderately or severely hemolyzed specimens are unsatisfactory because enzymes and intermediates liberated from the erythrocytes may affect the lag phase and the side reactions occurring in the assay system.

**EXPECTED VALUES**

Men	24 - 204 U/I	(0.39 - 3.40 µkat/l)
Women	24 - 173 U/I	(0.39 - 2.90 µkat/l)

Each laboratory should establish appropriate reference intervals related to its population.

**QUALITY CONTROL AND CALIBRATION**  
 It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:  
**QUANTINORM CHEMA**  
 with normal or close to normal control values  
**QUANTIPATH CHEMA**  
 with pathological control values.  
 If required, a multiparametric, human based calibrator is available:  
**AUTOCAL H**

Please contact Customer Care for further informations.

**TEST PERFORMANCE**  
**Linearity**  
 the method is linear up to 2000 U/I.  
 If a ΔA/min of 0.250 is exceeded, it is suggested to dilute sample 1+9 with saline solution and to repeat the test, multiplying the result by 10.

**Sensitivity/limit of detection (LOD)**  
 the limit of detection is 1.6 U/I.

**Interferences**  
 no interference was observed by the presence of:  
 hemoglobin ≤ 400 mg/dl  
 bilirubin ≤ 40 mg/dl  
 lipids ≤ 660 mg/dl

**Precision**

intra-assay (n=10)	mean (U/I)	SD (U/I)	CV%
sample 1	148.21	0.94	0.64
sample 2	464.75	3.98	0.86

inter-assay (n=20)	mean (U/I)	SD (U/I)	CV%
sample 1	148.35	1.33	0.90
sample 2	461.34	4.62	1.00

**Methods comparison**  
 a comparison between Chema CK-NAC FL and a commercially available product gave the following results:

CK NAC Chema = x  
 CK-NAC competitor = y  
 n = 100  
 $y = 1.04x - 3.10$   $r^2 = 0.9985$

**WASTE DISPOSAL**  
 This product is made to be used in professional laboratories.  
 P501: Dispose of contents according to national/international regulations.