APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ HITACHI 911/912	
TEST:	CHE
APP. CODE:	338
WAVELENGTH (Sec/Pri):	700 - 415
ASSAY:	RATE-A <i>TIME: 10</i> <i>POINT: 23-30</i>
SAMPLE VOL:	NORMAL: 5 DECREASE: 3 INCREASE: 8
	R1 VOLUME: 250 R2 VOLUME: 0 R3 VOLUME: 50 R4 VOLUME: 0
ABS LIMIT:	6500 - DEC
PROZONE LIMIT:	0 - UPPER
CALIB METHOD:	LINEAR (POINT: 2 - SPAN: 2 - WEIGHT: 0)
SD LIMIT:	0.250
DUPLICATE LIMIT:	3%
ST. 1 CONC:	0.0
EXPECTED VALUE:	4000 - 12000
UNIT:	U/I
INSTR. FACTOR (y=ax+b):	a=1 b=0

APPLICAZIONE / APPLICATION / APPLICATION / APLICACIÓN / ΠΡΟΓΡΑΜΜΑ OLYMPUS AU 400/480/600/640/680/2700 (Test code 868)	
TEST NAME:	CHE
SAMPLE:	Volume 5 µl Dilution 0 µl
REAGENTS:	R1 Volume 250 µl Dilution 0 µl R2 Volume 50 µl Dilution 0 µl
WAVELENGTH:	Pri. 410 Sec. 700
METHOD:	RATE
REACTION SLOPE:	-
MEASURING POINT 1:	First 17 Last 23
MEASURING POINT 2:	First Last
REAGENT OD LIMIT:	First L 1.0 First H 2.0 Last L 1.0 Last H 2.0
DYNAMIC RANGE:	L 430 H 25000
CORRELATION FACTOR:	A 1 B 0
LINEARITY LIMIT:	15%
UNIT:	U/I
CALIBRATION TYPE:	AB
FORMULA:	Y = AX + B

 Chema Diagnostica
Via Campania 2/4
60030 Monsano (AN) - ITALY - EU
phone +39 0731 605064
fax +39 0731 605672
e-mail: mail@chema.com
website: http://www.chema.com

ITALIANO rev. 07/07/2020

COLINESTERASI DGKC FL	
CH 2H120	5 x 20 + 2 x 10 ml
CH 6U120	2 x 50 + 2 x 10 ml

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INTENDED USE
Reagent for quantitative in vitro determination of cholinesterase in biological fluids.

PRINCIPLE OF THE METHOD
This reagent is formulated according to DGKC recommendations. Cholinesterase (pseudocholinesterase EC 3.1.1.8) catalyzes the hydrolysis of butyrylthiocholine, forming butyrate and thiocholine, which reduces the ferricyanide ions to ferrocyanide.
The decrease in absorbance is followed at 405 nm and it is proportional to cholinesterase activity in examined sample.

KIT COMPONENTS
For in vitro diagnostic use only.
The components of the kit are stable until expiration date on the label.
Keep away from direct light sources.

CHE R1 2H120 5 x 20 ml (liquid) white cap
6U120 2 x 50 ml (liquid) white cap

CHE R2 2H120 2 x 10 ml (liquid) red cap
6U120 2 x 10 ml (liquid) red cap

Composition in the test: sodium pyrophosphate 75 mM, pH 7.60, potassium ferricyanide 2 mM, butyrylthiocoline 15 mM, stabilizers.

Store all components at 2-8°C.

REAGENT PREPARATION
Use separate reagent ready to use.
Stability: up to expiration date on labels at 2-8°C.
Stability since first opening of vials: preferably within 60 days at 2-8°C -away from light sources-.
Caution: keep well refrigerated.

PRECAUTIONS
CHE R1: Danger. Causes serious eye damage (H318). Wear protective gloves. Eye protection (P280). IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing (P305+P351+P338). Immediately call a doctor (P310).

CHE R2: It is not classified as hazardous.

SPECIMEN
Serum, plasma (EDTA, heparinate only). Avoid hemolysis. ChE is stable in sample for at least 14 days whether the sample is stored at room temperature or under refrigeration.

EXPECTED VALUES
Total SChE:
Men: 5600 - 11200 U/I
Women: 4200 - 10800 U/I

Dibucaine number:
Normal homozygotes: > 75%
Heterozygotes: 35 - 75%
Atypical homozygotes: < 35%

Each laboratory should establish appropriate reference intervals related to its population.

QUALITY CONTROL AND CALIBRATION
It is suggested to perform an internal quality control. For this purpose the following human based control sera are available:

QUANTINORM CHEMA
with normal or close to normal control values
QUANTIPATH CHEMA
with pathological control values.
If required, a multiparametric, human based calibrator is available:

AUTOCAL H
Please contact Customer Care for further information.

TEST PERFORMANCE
Linearity
the method is linear up to 25000 U/I.
If a ΔA/min of 0.30 is exceeded, it is suggested to dilute sample 1+9 with saline solution and to repeat the test, multiplying the result by 10.

Sensitivity/limit of detection (LOD)
the limit of detection is 432.3 U/I.

Interferences
no interference was observed by the presence of:
hemoglobin ≤ 500 mg/dl
bilirubin ≤ 40 mg/dl
lipids ≤ 800 mg/dl

Precision
intra-assay (n=10) mean (U/I) SD (U/I) CV%
sample 1 5972.9 122.8 2.1
sample 2 5743.8 57.5 1.0

inter-assay (n=20) mean (U/I) SD (U/I) CV%
sample 1 5808.4 113.4 2.0
sample 2 5753.5 99.6 1.7

Methods comparison
a comparison between Chema Cholinesterase FL and a commercially available product gave the following results:

SChE Chema = x
SChE competitor = y
n = 107

y = 0.985x + 51.7 U/I r²=0.996

WASTE DISPOSAL
This product is made to be used in professional laboratories.
P501: Dispose of contents according to national/international regulations.

CHOLINESTÉRASE DGKC FL	
CH 2H120	5 x 20 + 2 x 10 ml
CH 6U120	2 x 50 + 2 x 10 ml

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